


SATURN FOLDING WALKER



SATURN FOLDING WALKER - WHEELS



BASIK AND INVACARE RECIPROCAL / AVENCIA WALKER



NAME	SATURN FOLDING WALKER	SATURN FOLDING WALKER – WHEELS	BASIK AND INVACARE RECIPROCAL / AVENCIA WALKER
PRODUCT CODE	WAL001	WAL002	P402/ALK713L
FRAME	ALUMINIUM RAW FINISH	ALUMINIUM RAW FINISH	
FERRELL	22MM	22MM	
ADJUSTABILITY	SINGLE PIN, FOLDING INWARDS	SINGLE PIN, FOLDING INWARDS	
GRIP	STD FOAM GRIP	STD FOAM GRIP	
SIZE	ADULT	ADULT	78-95 CM HEIGHT
GROSS WEIGHT	2.5KG	2.5KG	6 KG
COLOUR	GREY	GREY	
USER WEIGHT	MAX 100KG GROSS WEIGHT – 3KG	MAX 100KG GROSS WEIGHT – 3KG	NETT WEIGHT 2.3 KG

INVACARE WALKLITE



Flip-up seat allows user to step inside the frame



Model no. 6300-JRA (Junior to Adult)

Model no. 6300-ATA (Adult to Tall)

Folds into a compact unit for storage. Built-in seat handle allows for easy transport

PRODUCT CODE: 6300ATA

Features

- Lightweight aluminium frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanisms
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Rear spring-loaded Invacare® SureGlide™ brakes
- Adjustable seat-to-floor height

INVACARE HEMI SIDESTEPPER & FORIA TWO STEP



PRODUCT CODE: 6252

Invacare® Adult Hemi Walker / Sidestepper Model no 6252-Hemi Walker

- Wider base provides more stability than a quad cane
- 14.75" x 17.75" width at high height setting
- 12.25" x 16.5" width at low height setting



PRODUCT CODE: P421B

New Invacare® Two Step Walker Model no P421B

- Combination single release folding walker and sit to stand assistance device
- Can also be used as a toilet safety frame
- Adjustable height settings
- Accessories currently not available

4 WHEEL ALUMINIUM CABLE BRAKE ROLLATOR



INVACARE BANJO ROLLATOR



NAME	4 WHEEL WALKER WITH CABLE BRAKES	INVACARE BANJO ROLLATOR
PRODUCT CODE	10279	1452442
FRAME	ALUMINIUM	MILD STEEL
CASTOR	FOUR 6"X1" SOLID CASTORS	8 X 2" SOLID
BRAKE	CABLE BRAKE	CABLE BRAKE
STD ACCESSORIES	BASKET, VINYL UPHOLSTERED SEAT	BASKET, CANE HOLDER, PERSPEX TRAY
FOLDING	YES	YES
HANDLES	ERGONOMIC GRIP & HEIGHT ADJUSTABLE ARM	ERGONOMIC GRIP & HEIGHT ADJUSTABLE
WEIGHT CAPACITY	LOAD CAPACITY – 100KG GROSS WEIGHT – 8KG (EACH)	LOAD CAPACITY – 130KG GROSS WEIGHT – 9KG
COLOUR	BLUE, MAROON	SILVER
PACKING	38X60X87CM (2 PER BOX)	45X60X87CM

INVACARE DOLOMITE JAZZ 610 ROLLATOR



reddot design award winner 2013



Big basket
Provides extra storing space and can be closed off. Available in grey or red.



Back strap
Comfortable and safe support to lean against while sitting. Available with (top picture) or without extra padding (lower picture).



Cane holder
Allows to safely transport a cane or a crutch and makes it easily accessible.

INVACARE ROLLITE



Folds into a compact unit for storage. Build-in seat handle allows for easy transport

Colour: Midnight Blue
Model No's:
68100 Rollite Adult w/8" wheels,
68100-TA Rollite Tall Adult w/8" wheels (shown)

Features

- Lightweight aluminium frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanism
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Ergonomic handbrakes lock the rear wheels for security
- Adjustable seat-to-floor height

Linking style with function

The DOLOMITE JAZZ is an innovative lightweight rollator designed to offer excellent support and stability, both during indoor and outdoor walks. This new stylish rollator includes numerous smart features making it the ideal walking companion. The latest member of the DOLOMITE family combines style with function. Its ergonomic design allows a natural walking position inside the rollator.

NAME	INVACARE DOLOMITE JAZZ 610 ROLLATOR	68100-TA ROLLITE TALL ADULT W/8" WHEELS (MIDNIGHT BLUE)
PRODUCT CODE	JAZZ610	68100TA
DEPTH OPENED AT BASE	690MM	26"
DEPTH FOLDED		7"
WIDTH INSIDE HANDGRIPS		18.25"
WIDTH INSIDE BASE OF LEGS	595MM	22"
WIDTH OVERALL	690MM	27"
HEIGHT ADJUSTMENT MIN – MAX	750-1000MM	35-39"
*APPROX CLIENT HEIGHT MIN – MAX		5'11"-6'6"
SEAT-TO-FLOOR HEIGHT MIN – MAX	610MM	22"-26"
SEAT WIDTH X DEPTH		17"X9.5"
WHEELS	200MM	8"
PRODUCT WEIGHT	7.8KG	16.5LB
WEIGHT CAPACITY		300LB
UNITS / CARTON	1	1
LIMITED WARRANTY		
BRAKES		5 YEARS
OTHER COMPONENTS**		1 YEAR

ELBOW CRUTCH



INVACARE CRUTCH



Invacare® Advantik Advantik is a range of stable light weight crutches with a weight capacity of 130kg. The click-together function is convenient when quickly needing to get one hand free. Reflector is standard and XL.



Invacare® Basik


Load Capacity on Basik of 100Kg

NAME	ALUMINIUM ELBOW CRUTCH	INVACARE CRUTCH
PRODUCT CODE	CRU001	CRU003
FRAME	ALUMINIUM	
FERRELL	22MM	
ADJUSTABILITY	SINGLE PIN	
GRIP	STD CONTOURED GRIP	
SIZE	KIDDIES, YOUTH, ADULT	
GROSS WEIGHT	800G	
COLOUR	GREY	
USER WEIGHT	MAX 100 KG	MAX 100 KG

CANES & STICKS



FOLDING CANES



UNDERARM CRUTCH



NAME	QUAD CANE /TRIPOD CANE	ADJUSTABLE FOLDING WALKING STICK	ALUMINIUM UNDERARM CRUTCH
PRODUCT CODE	ALK521	APC20061	7021(M) 7022(L)
FRAME	ALUMINIUM & MILD STEEL	ALUMINIUM	ALUMINIUM
FERRULE	15MM	19MM	22M
ADJUSTABILITY	SINGLE PIN 75 – 95 CM	SINGLE PIN IN 1" INCREMENTS FOLDED – 30CM ADJUSTABLE FROM 85CM – 96CM	PATIENT HEIGHT ADJ HEIGHT 158-178(CM)112-132(CM) 178-198(CM) 132-152 (CM)
GRIP	STD FOAM GRIP	STD CONTOURED GRIP	FOAM
SIZE	ADULT	ADULT	
GROSS WEIGHT	1.8KG, PACKING: 26X24X78CM	400G	LOAD CAPACITY – 100KG GROSS WEIGHT – 1KG EA
COLOUR	BLACK	BRONZE – VARIOUS PATTERNS & COLOURS (STOCK DEPENDENT)	GREY
USER WEIGHT	MAX 110KG	MAX 100KG	MAX 100KG (120KG - TALL)
BASE FOR QUAD CANE	LARGE: 33X40X75CM		

REHABILITATION WALKER



The newest model of walking supports
Easy-Walker

Connecting people with Health,
Happiness and Dignity ...



PRODUCT CODE: EAS001