

**BARIATRIC ROLLATOR**



**Features**

- Flip-up 2" foam padded seat w/ zipper compartment
- Removable padded backrest
- Folds for ease of storage and transport
- Height adjustable loop handbrakes with ergonomic handles
- Indoor/Outdoor 7.5" wheels

**HEAVY DUTY PADDLE WALKER**



**Features**

- 500lb weight capacity
- Extended frame for extra stability
- Long-lasting PVC handgrips

NAME	BARIATRIC ROLLATOR	HEAVY DUTY PADDLE WALKER
PRODUCT CODE	68100	6291-HDA
MODEL NO	66550	6291-HDA
DEPTH OPENED AT BASE	29.5"	17.5"
WIDTH INSIDE HANDGRIPS	22"	20"
WIDTH INSIDE BASE OF LEGS	20.25"	24"
WIDTH OVERALL	27"	
HEIGHT ADJUSTMENT (MIN – MAX)	33.75"-38"	
SEAT WIDTH X DEPTH	18"x14"	
PRODUCT WEIGHT	26lb	6lb.11oz
WEIGHT CAPACITY	500lb	500lb
DEPTH FOLDED		4"
HEIGHT ADJUSTMENT (MIN-MAX) (ADJUSTS IN 1° INCREMENTS)		30.38"-37.38"
APPROXIMATE PATIENT HEIGHT (MIN-MAX)		5'0"-6'3"
QUANTITY		2 per carton
LIMITED WARRANTY		5 years

INVACARE BARIATRIC CANE



INVACARE BIG50 CRUTCH



Features

- 700lb weight capacity
- Right- or left hand use
- Height adjustable
- Comfortable PVC hand grip
- Locking ring for additional stability
- Durable steel construction

Invacare® Big 150

- A sturdy crutch with weight capacity for users up to 150kg.
- The 1080mm crutch can be cut to fit any user.
- Adjustable arm support.

NAME	INVACARE QUAD CANE	INVACARE BARIATRIC CANE	INVACARE BARIATRIC CANE
PRODUCT CODE	7835-2	8940-6	152986-8
MODEL NO	7831-2	7835-2	8940-6
BASE	SMALL	LARGE	OFFSET
COLOUR	BLACK	BLACK	BLACK
ADJUSTMENT (MIN – MAX)	29.50"-38.50"	29.50"-38.50"	29"-38"
CLIENT HEIGHT (APPROX)	4'11"-6'5"	4'11"-6'5"	4'10"-6'4"
FOOTPRINT	12"-8.50"	12"X8.50"	-
WEIGHT	3LB	4LB	3LB
WEIGHT CAPACITY	700LB	700LB	700LB
QUANTITY	2 PER CARTON	2 PER CARTON	6 PER CARTON
LIMITED WARRANTY	3 YEARS	3 YEARS	3 YEARS

**BARIATRIC SHOWER CHAIR**



**BARIATRIC TRANSFER BENCH**



**Features**

- 700 lb weight capacity
- Wide, deep seating surface
- Durable aluminium frame with dual cross brace for added stability
- Backrest easily snaps into assembled seat frame with no tools required
- Contoured seat and backrest
- Adjustable height legs

NAME	BARIATRIC SHOWER CHAIR	BARIATRIC TRANSFER BENCH
PRODUCT CODE	9785	9675
MODEL NO	9785-2	9675-2
SEAT HEIGHT (MIN – MAX)	16.25"-20"	16.25" – 20.25"
SEAT DEPTH	16.25"	16"
SEAT WIDTH	17"	17.25"
SEAT WIDTH (INCLUDING TRANSFER SECTION)		27"
SEAT WIDTH (INSIDE TUB WITH HANDRAIL)		19"
BACK HEIGHT FROM SEAT	13.50"	
OVERALL WIDTH AT BASE (MIN – MAX)	18.25" – 18.75"	
OVERALL DEPTH AT BASE (MIN – MAX)	21" – 22"	
OVERALL HEIGHT		28" – 32.25"
PRODUCT WEIGHT	11LB	16LB
WEIGHT CAPACITY	700LB	700LB
QUANTITY	2 PER CARTON	2 PER CARTON
LIMITED WARRANTY	3 YEARS	3 YEARS

**I-CLASS HEAVY DUTY TRANSFER BENCH**



**Features**

- 400lb weight capacity
- Study, durable, one-piece blow-moulded seat
- Textured seat with drain holes
- Ergonomically shaped reversible backrest snaps in and out easily without tools
- Built-in soap dish and hand-held shower holder

**BARIATRIC DROP-ARM COMMODE**



**Features**

- 1000lb weight capacity
- Easy-to-operate drop-arm mechanism
- Extra-wide snap-on seat
- Removable padded backrest
- Adjustable height legs
- Durable steel frame design
- Comes complete with pail and lid

NAME	i-CLASS HEAVY DUTY TRANSFER BENCH	NAME	BARIATRIC DROP-ARM COMMODE
PRODUCT CODE	96704	PRODUCT CODE	6599
MODEL NO.	9670, 9670U, 9670C	MODEL NO.	6599
SEAT HEIGHT (MIN – MAX)	18:-22.5"	BASE WIDTH OVERALL	28"
SEAT DEPTH	18"	WIDTH BETWEEN ARMS	25"
SEAT WIDTH	33.25"	BASE DEPTH OVERALL	20.5"
OVERALL HEIGHT MIN – MAX	31.25:-35.75"	SEAT SURFACE DEPTH	18"
CLEARANCE BENEATH FRAME	14.5"-31"	SEAT SURFACE WIDTH	23.5"
OVERALL WIDTH AT BASE	30.5"-31"	SEAT OPENING WIDTH	8.25"
OVERALL DEPTH AT BASE	16.25"-17"	SEAT OPENING DEPTH	10.25"
HEIGHT (ADJUSTMENTS IN .5° INCREMENTS)	9	HEIGHT OVERALL	27"-31"
PRODUCT WEIGHT	12lb	SEAT-TO-FLOOR HEIGHT	19.25"-23.25"
WEIGHT CAPACITY	400lb	SEAT-TO-ARM HEIGHT	7.75"
QUANTITY PER CARTON		PRODUCT WEIGHT	21lb
9670	2	WEIGHT CAPACITY	1000lb
9670U	1	LIMITED WARRANTY	5 years
9670C	2		
LIMITED WARRANTY	3 years		

**BARIATRIC PATIENT HOIST ELECTRIC**



**BARIATRIC TRAPEZE**



**Features**

- 1000 lb client weight capacity
- Three locking overhead boom angle positions for comfortable, secure use
- Tool-free design simplifies assembly
- Four piece sub-assembly significantly reduced setup time
- Lightweight design for easy delivery
- Free-standing structure permits use with Bariatric, Home Care or non-medical bed

<b>NAME</b>	ELECTRIC-BARIATRIC PATIENT HOIST
<b>PRODUCT CODE</b>	SMA10601E
<b>SLING</b>	BARIATRIC
<b>FRAME</b>	POWDER COATED MILD STEEL
<b>BASE</b>	66CM – 115CM U BASE – ADJUSTABLE
<b>CASTOR</b>	FOUR 4" SOLD CASTORS – REAR WITH BRAKES
<b>CONTROL</b>	LINAK ACTUATOR
	OVERRIDE FOR MANUAL RELEASE IN POWER OUT
<b>CARRY CAPACITY</b>	WEIGHT – 215KG GROSS WEIGHT – 65KG
<b>HEIGHT ADJUSTMENT</b>	865 – 2075MM
<b>CHARGE TIME</b>	4 TO 12 HOURS

<b>NAME</b>	BARIATRIC TRAPEZE
<b>PRODUCT CODE</b>	BARTRAP
<b>MODEL NO.</b>	BARTRAP
<b>OVERALL HEIGHT</b>	73"
<b>REACH</b>	35.5"
<b>BASE WIDTH</b>	40"
<b>BASE LENGTH</b>	50.5"
<b>PRODUCT WEIGHT</b>	72lb
<b>SHIPPING WEIGHT</b>	80lb
<b>WEIGHT CAPACITY</b>	1000lb
<b>LIMITED WARRANTY</b>	3 years
<b>HCPSC CODE</b>	E0912