

INVACARE ORCA AND ORCA XL BATHLIFT



Lift in operation

- Maintenance free actuator, powerful and quiet
- Removable and machine-washable covers
- High backrest with contoured neckrest, can be positioned and folded with just one hand
- Optimised base plate – significantly more space in the bathtub
- Patented, self-released suction cups: can be easily removed
- Backrest can be reclined up to 40°
- Wide range of accessories for individual customisation
- The ergonomic, floating hand control fits securely in the hand and is easy to grasp



Can also be used without cover



Can be disassembled into 2 parts



High-performance actuator

PRODUCT CODE ORCA: 1470322
 PRODUCT CODE ORCA XL: 1471066

Technical Data

Order Number	Aquatec Orca	Aquatec Orca XL
Bathtub lift with cover	1470322	1471066
Bathtub lift without cover	1470323	
Dimensions		
Seat width without side flaps	37.5 cm	37.5cm
Seat width with side flaps	71.5cm* (A)	71.5 cm (A)
Seat depth	49 cm (B)	49 cm (B)
Seat height max	42 cm (C)	42 cm (C)
Seat height min	6 cm	6cm
Width of backrest	36 cm (D)	36 cm (D)
Height of backrest	66.5 cm (E)	66.5 cm (E)
Total depth (backrest upright)	65 cm (F)	65 cm (F)
Total depth (backrest reclined)	92 cm	92 cm
Total height with backrest	107.5 cm (G)	107.5 cm (G)
Base area	58 x 28.5 cm	58 x 28.5 cm
Cover thickness	0.3 cm	0.3 cm
More Information		
Total weight without hand control	121.1 kg	12.1 kg
Total weight without battery pack	-	-
Weight, seat	7.5 kg	7.5 kg
Weight, backrest	4.6 kg	4.6 lg
Weight, hand control	0.4 kg	0.4 kg
Weight, battery pack	-	-
Load capacity	135 kg	170 kg
Backrest can be reclined by up to	40°	40°
Colour	White	White

MANGAR BATH CUSHION



Benefits / Features

- Exceptionally stable
- Lifts up to 150kg (24 stone)
- Simple to use
- Quick and easy to clean
- Battery powered – no dangerous power cables
- Portable with soft carry handle
- Lightweight – weight only 2kgs (4½ lbs)
- Easy to remove from the bath and compact for storage
- Comfortable for frail and underweight individuals
- Battery level indicator light and built-in safety feature preventing user from being lowered into the bath if there is not enough battery power to raise them back up again
- Minimal number of seams to protect fragile skin
- Can be kept slight inflated to provide lumbar support at the bottom of the bath
- Powered by our Airflo compressor – only one compressor required to multiple pneumatic products

PRODUCT CODE: MANGAR

The New Bathing Cushion is designed to offer increased stability and comfort. Weighing only 2kgs (4.4lbs) and featuring a new carry handle, the Bathing Cushion is completely portable. It can be set up or removed in seconds and requires no fixtures to an existing bathroom, making it ideal for taking away on visits or on holiday.

- Our most stable inflatable bathlift ever
- Now even easier to transport
- Will not lower you into the bath unless there is sufficient battery power to raise you back out again
- Compact and easy to store away
- Lightweight
- Supportive backrest
- Wide seat for comfort and stability
- Fewer seams to protect fragile skin



ADJUSTABLE SHOWER CHAIR



ADJUSTABLE SHOWER CHAIR - WITH TRANSFER



Q/R ADJUSTABLE SHOWER CHAIR - WITH ARMS

Model	Size (mm)	48	75	40-50	74-84	40	36	70x52x47	N.W (kg)	PCS CN	G.W (kg)
ALP540L									4	2	10



SHOWER CHAIR



ALK405L-3" Shower Chair

FOLDING ALUMINIUM SHOWER CHAIR

Model	Size (mm)	64	41	56	91	38	36	109x56x32	N.W (kg)	PCS CN	G.W (kg)
ALK405L-3"									7	2	16



NAME	ADJUSTABLE SHOWER CHAIR – WITH & WITHOUT BACK	Q/R ADJUSTABLE SHOWER CHAIR W/ARM	NAME	FOLDING ALUMINIUM SHOWER CHAIR
PRODUCT CODE	APC5001 & 10560	5041 & 5042	PRODUCT CODE	5017WADJ
SEAT	480X300 (MM)	520X400 (MM)	DIMENSION	50X54X90-98CM
ADJUSTABLE	420X520 (MM)	705-805 (MM) & 580-680 (MM)	ADJUSTABLE	4 LEVEL ADJUSTABLE
WEIGHT	2KG (NO BACK) 3 KG (WITH BACK)	2KG (NO BACK) 3KG (WITH BACK)		
LOAD CAPACITY	90 KG	105KG		
OPTIONS	ALSO AVAILABLE WITH "U" SHAPE SEAT SEE INSERT PICTURE		MATERIAL	1" ALUMINIUM FRAME, PE BOARD, TPR TIP, POWDER COATED

BASIK FOLDING SHOWER SEAT



FOLDING SHOWER SEAT



FOLDING SHOWER SEAT - WITH LEGS



DELUXE SHOWER SEAT



NAME	BASIK FOLDING SHOWER SEAT	FOLDING SHOWER MOULDED SEAT	FOLDING SHOWER SEAT – WITH LEGS	DELUXE FOLDING SHOWER SEAT
PRODUCT CODE	SHM5054	SMN6104	SMN5051	FSC
FRAME	POWDER COATED STEEL	POWDER COATED MILD STEEL	POWDER COATED MILD STEEL	STAINLESS STEEL
SEAT SIZE	37X23X2.7CM	490X430X480 (MM)	400X400 (MM)	460X390 (MM)
SEAT	BLOW MOULDED PLASTIC SEAT	ONE PIECE MOULDED PLASTIC SEAT	ONE PIECE MOULDED PLASTIC SEAT	ONCE PIECE SOLD PNENOLIC SEAT
WEIGHT		3.6KG, GROSS WEIGHT 4.2KG	4KG, GROSS WEIGHT 4.5KG	7KG, GROSS WEIGHT 7.5KG
LOAD CAPACITY	300LBS	80KG	100KG	90KG
PACKING	38X33X6CM	SEAT 37DX40D (CM) 50CM FROM WALL TO FRONT OF SEAT	490X600X150 (MM)	580X450X150 (MM)

SHOWER STOOLS AND SHOWER CHAIRS



Aquatec Pico

Aquatec Pico 3 in 1

Technical Data


PRODUCT CODE: PICO & PICO3

Order Number	Aquatec Pico
Shower stool	9.20.001
Shower stool with armrests	9.21.001
Shower chair with backrest	9.23.001
Shower chair complete	9.22.001
Backrest	16469
Armrests (pair)	16468

Dimensions	
Seat width	43 cm (A)
Seat depth	42 cm (B)
Seat height	42.5-57.5 cm (C)
Total width	57.5 cm (D)
Total depth	52 cm (E)
Total height	82.5-97.5 cm (F)
Height to armrests	60.5-75.5 cm (G)
Width between the armrests	43 cm (H)
Width, hygiene recess	18 cm
Depth, hygiene recess	14.5 cm
Graduation, height adjustment	2.5 cm

More Information	
Weight, shower stool	3.2 kg
Weight, shower stool with armrests	4.1 kg
Weight, shower chair with backrest	4.0 kg
Weight, shower chair complete	4.8 kg
Load capacity	160 kg
Colour	White

Invacare I-Fit 9781E



Invacare I-Fit 9781E

- Height Adjustable from 38cm to 55 cm
- Removable backrest
- Ergonomic shaped seat with backrest for extra comfort
- Easy to clean
- Light weight and easy to transport



Technical Data

PRODUCT CODE: 9781e

Order Numbers	Invacare I-Fit 9781E
Shower Stool	
Shower Chair	9781E

Dimensions	
Seat width	34
Seat depth	37
Seat height	35-55
Total width	50
Total depth	43
Total height	76-93

More Information	
Weight	4.5 kg
Load Capacity	180 kg
Colour	White

INVACARE ALIZE AND SWIVEL BATHER WITH LOCK



INVACARE FOLDING MOBILE SHOWER CHAIRS



INVACARE® Folding Shower Chair With Back Features
 Model #9981
 Wide contoured seat and back
 Ergonomically positioned drain holes allow consistent drainage
 Hand-held show holders are moulded in for convenience
 Non-marring rubber tips

NAME	ALIZE & SWIVEL BATHER WIT LOCK	INVACARE FOLDING MOBILE SHOWER CHAIR
PRODUCT CODE	SMH-5025 & H3000	9981
SEAT	PLASTIC MOULDED SEAT WITH U CUT OUT 400X400 (MM) BETWEEN ARMS – 530 (MM)	
OPTIONS	SWIVEL – WITH LOCK	
WEIGHT	GROSS WEIGHT – 7.5KG	
LOAD CAPACITY	100KG & 130KG FOR ALIZE	
PACKING	720X540X240 (MM)	

SLATTED BATH BENCH



ADJUSTABLE BATH BOARD



NAME	SLATTED BATH BENCH	ADJUSTABLE BATH BOARD
PRODUCT CODE	BAT004	10440
SEAT	ADJUSTABLE SLATTED MOULDED PLASTIC	ADJUSTABLE ONE PIECE MOULDED PLASTIC
PACKING	480X340X220 (MM)	750X340X140 (MM)
SEAT SIZE	45X300 (MM) HEIGHT 180	700X310 (MM)
WEIGHT	2.3KG	2.8KG
LOAD CAPACITY	GROSS WEIGHT 3KG 90KG	GROSS WEIGHT 3.3KG 100KG
OPTIONS		1) WITH TRANSFER DISK 2) SLATTED BASE

INVACARE BATH MAT



PRODUCT CODE: H190



Model #	Overall Dimensions	Product Weight	Weight Capacity	Units Per Carton	Limited Warranty
H190	.4" H 16" W 32" D	3.75 lb.	286 lb.	6	1 Year

Features

Module # H190

- High-quality bathtub mat
- Held securely in the bathtub by numerous small suction cups on the underside of the mat
- Non-slip surface for a secure hold while bathing or showering
- Special contour design to fit around drain
- Latex free
- Extra long size

INVACARE MARINA



PRODUCT CODE: 1541330

- Attractively designed bathtub board
- Non-slip seat surface
- Hand grip provides additional safety
- Practical shower head holder on both sides
- Individually adjustable using adjustment elements
- Also available as an XL version
- Load Capacity 150kg

AQUATEC TRANSFER AIDS



Aquatec Step

PRODUCT CODE: 3.00.002

- Stable frame
- Feet with non-slip rubber caps
- Non-slip, removable cover, machine-washable

INVACARE H160 GRABHANDLES



- Grab handle for the bathtub
- Provides additional safety when getting into and out of the bath
- Easy assembly using quick release screws no tools required
- Safely secured on the edge of the bath using a clamping mechanism made of steel and a pressure plates with anti-slip surface