

STANDARD COMMODE



FOLDING HOLIDAY COMMODE



NAME	STANDARD STEEL COMMODOE	FOLDING HOLIDAY COMMODOE
PRODUCT CODE	Cm01	HCOM
SIZE	18"	18"
SEAT TO FLOOR	ADJUSTABLE 37 – 52 CM	ADJUSTABLE 50 – 60CM
FRAME	MILD STEEL CHROME / POWDER COATED	MILD STEEL POWDER COATED
COLOUR	CHROME / GREY	GREY
INCLUDES	PLASTIC TOILET SEAT REMOVABLE BUCKET LID FOR TOILET SEAT	PLASTIC TOILET SEAT REMOVABLE BUCKET LID FOR TOILET SEAT
LOAD CAPACITY	100KG	85KG
WEIGHT	WEIGHT – 6KG EACH GROSS WEIGHT – 7 KG	WEIGHT – 6KG EACH GROSS WEIGHT – 14KG (2 PER BOX
FOLDED		84X52X14CM (SEAT AND BUCKET REMOVED)

SATURN TAXI COMMODE



SHOWER S.P. COMMODE



NAME	SATURN TAXI COMMODE	SATURN SHOWER COMMODE
PRODUCT CODE	COMTAXFR	COMALU
SIZE	18"	18"
FRAME	CHROME PLATED STEEL	ALUMINIUM WHITE POWDER COATED
FOOTREST	REMOVABLE	REMOVABLE
CASTOR	5" RUBBER CASTOR – REAR WITH BRAKES	5X1" SOLD
REAR WHEEL FOR SP	24" MAG WITH 24X1" SOLID TYRE	
REAR RIMS		24X13/8 MAG RIM WITH 24X13/8 PU GREY TYRE
UPHOLSTERY	VINYL BACK	VINYL
CARRY CAPACITY	WEIGHT – 115KG GROSS WEIGHT – 19KG	WEIGHT – 115KG GROSS WEIGHT – 17KG
SEAT	POLYURETHANE – 18"	POLYURETHANE – 18"
ARMS	FLIP DOWN FULL ARMS – TAXI FLIP UP FULL ARMS – SP	REMOVABLE & HEIGHT ADJUSTABLE T ARMS
BACK	REMOVABLE FOR TRANSIT	REMOVABLE FOR TRANSIT
PACKING	TAXI – 600WX600DX550H (MM)	600X610X630 (MM)
SEAT-TO-FLOOR HEIGHT	SP – 700WX830DX660H (MM) SP – 19" TAXI – 20"	

FOLDING DELUXE COMMODE



FANTASTIC COMMODE FOR TRAVELING AND OUTDOOR ACTIVITY!!!



CARRY BAG SIZE HANDY



FOLDING



HANDY TIGHTEN KNOB

NAME	FOLDING DELUXE COMMODE
PRODUCT CODE	SMA-LFC
OPERATION SIZE	57 X 53 X 70-83 CM
FOLDING SIZE	62 X 40 X 23 CM
BAG SIZE	67 X 43 X 26 CM

BASIK COMMODES



ALK613

Model	Size	Icon 1	Icon 2	Icon 3	Icon 4	Icon 5	Icon 6	Icon 7
ALK613	56	-	54	75	75	45-55	72-82	50


44	35	100	50x22x72	8	10
----	----	-----	----------	---	----




ALK699L

Model	Size	Icon 1	Icon 2	Icon 3	Icon 4	Icon 5	Icon 6	Icon 7
ALK699L	54	-	54	10	10	54	80	55

44	35	100	80x70x62	7	10
----	----	-----	----------	---	----



ALK612



ALK615

Model	Size	Icon 1	Icon 2	Icon 3	Icon 4	Icon 5	Icon 6	N.W	PCS	G.W
ALK612	50	56	45-55	72-82	44	35	50x22x72	7	1	9
ALK615	48	44	48	82	43	34	50x25x82	9	1	11



ALK694

Model	Size	Icon 1	Icon 2	Icon 3	Icon 4	Icon 5
ALK694	55	-	43	12.5	12.5	

54	95	93	43	37	100	65x65x66	M	N
----	----	----	----	----	-----	----------	---	---

INVACARE H440 OMEGA ECO



PRODUCT CODE: 1542903

Technical Data

Order Numbers	H440 Omega Eco
Commode Chair, 5 pack	1542903
Dimensions	
Seat width	450 mm
Seat depth	430 mm
Seat height	405 – 585 mm
Total width	580 mm
Total depth	530 mm
Total height	790 – 970 mm
Height to armrests	620 – 800 mm
Height, ground to lower seat edge:	
With holder for toilet pan	335 – 515 mm
Without holder for toilet pan	345 – 525 mm
Height floor to seat plate	405 585 mm
Width between the armrests	470 mm
Width, cut out in seat plate	200 mm
Depth, cut out in seat plate	255 mm
More Information	
Load capacity	140 kg
Weight (kg)	7.5 kg
Colour	Lavender blue

INVACARE H750 COMMODE



Invacare Omega H750

- Compact modern design
- Castors with brakes for safety
- Padded back, seat and arms
- Integral foot support

PRODUCT CODE: 1471044

Order Numbers	Invacare Omega H750
Commode chair, mobile	
Commode chair, mobile, in special construction	
Commode chair, mobile, lavender blue	1471043
Commode chair, mobile, blue	1471044
Commode chair, lavender blue	
5 pack	
Commode chair, blue	
5 pack	
Commode chair, lavender blue, 2 pack	
Dimensions	
Seat width	50 cm (A)
Seat depth	46 cm (B)
Seat height	48 cm (C)
Total width	63.5 cm (D)
Total depth	79 cm (E)
Total height	92 cm (F)
Depth without footrests	
Footplate depth	
Height, footrests to seat	
Height, ground to lower seat edge:	
with holder for toilet pan	36 cm (K)
without holder for toilet pan	-
Height, floor to cover plate	50.5 cm (L)
Width between the armrests	51 cm (M)
Width, cut out in seat plate	21.5 cm
Depth, cut out in seat plate	25.5 cm
Ø Castors	100 mm
More Information	
Weight	15.0 kg
Load capacity	130 kg
Colour	Lavender blue / Blue

INVACARE LIMA COMMODE



Lima with 24" wheels

PRODUCT CODE: H273

INVACARE AQUATEC OCEAN VIP



AQUATEC



PRODUCT CODE: 1470712

Technical data

H263	58 cm	47 cm	57 cm	104,5 cm
H273	58 cm	47 cm	68 cm	111,5 cm
H263	100 cm	21 cm	14 kg	135 kg
H273	100 cm	30 cm	18,5 kg	135 kg



Seat tilt
Tilt angle 0° - 35°.



Headrest
Height and depth-adjustable headrests included.



Adjustable seat height
Adjustable in 3 steps, no tools required.



Directional lock
Easy to operate brakes on 3 castors, 1 castor with directional lock to help push the chairs down straight floors.



Wider XL armrests
Increases the distance between the armrests by 4 cm on each side.
Retrofit kit for *Ocean Vip*, *Ocean Dual Vip* and *Ocean E-Vip*
1470086 Wider XL armrests (pair)



Ergonomic special soft seat
Extra soft, textured surface
1531259 soft for *Ocean*, *Ocean XL*, *Ocean Vip*, *Ocean Vip XL*, *Ocean Dual Vip*