

SATURN COASTER



SATURN HERCULES



USE: Active Med-Active Passive
 ENVIRONMENT: Indoor Indoor/Outdoor Outdoor/Indoor
 ACTIVITY LEVEL: High Medium Low
 ADJUSTABILITY: Folding Back-height Back angle
 Seat angle COG Q/R wheel
 CONFIGURABILITY: TAB Full arm E/leg rest
 H/Adj arm 8 x 2" S/depth

USE: Active Med-Active Passive
 ENVIRONMENT: Indoor Indoor/Outdoor Outdoor/Indoor
 ACTIVITY LEVEL: High Medium Low
 ADJUSTABILITY: Folding Back-height Back angle
 Seat angle COG Q/R wheel
 CONFIGURABILITY: TAB Full arm E/leg rest
 H/Adj arm 8 x 2" S/depth

| NAME | SATURN COASTER | SATURN HERCULES |
|----------------------|--|--|
| PRODUCT CODE | COA | HER |
| SIZE | 14", 16", 18", 20" | 14", 16", 18", 20" |
| COLOURS | COBOLT BLUE | COBOLT BLUE |
| FRAME | REINFORCED POWDER COATED MILD STEEL | REINFORCED POWDER COATED MILD STEEL |
| ARMS | REMOVABLE DESK | REMOVABLE DESK OR FULL ARM |
| FOOTREST | REMOVABLE, ELEVATING LEG REST UPCHARGE OPTION | REMOVABLE, ELEVATING LEG REST OPTIONAL EXTRA |
| CASTORS | 8x1" SOLID CASTOR ON STEEL FORK | 8x2" PU CASTOR ON 5MM STEEL FORK |
| RIMS | 24X1 ^{3/8} MAG RIM WITH 24X1 ^{3/8} PU TYRE - FIXED AXLE | 24X1 ^{3/8} MAG RIM WITH 24X1.75 PU TYRE - FIXED AXLE |
| UPHOLSTERY | BLACK REINFORCED NYLON | BLACK REINFORCED NYLON |
| STRENGTH | SABS COMPLIANT AND TESTED, CHROME PUSHRIMS, STEEL SIDE PANELS, DOUBLE CROSSBARS ON THE 20" MODEL | SABS COMPLIANT AND TESTED, CHROME PUSHRIMS, STEEL SIDE PANELS, DOUBLE CROSSBARS ON THE 20" MODEL |
| LOAD CAPACITY/WEIGHT | LOAD CAPACITY 120KG WEIGHT - 20KG | LOAD CAPACITY 125KG WEIGHT - 21KG |
| PACKING | 100 X 35 X 85 (CM) | 100 X 35 X 85 (CM) |

SATURN PORTERS



OMEGA ECONOMY



- USE: Active Med-Active Passive
 ENVIRONMENT: Indoor Indoor/Outdoor Outdoor/Indoor
 ACTIVITY LEVEL: High Medium Low
 ADJUSTABILITY: Folding Back-height Back angle
 Seat angle COG Q/R wheel
 CONFIGURABILITY: TAB Full arm E/leg rest
 H/Adj arm 8 x 2" S/depth

- USE: Active Med-Active Passive
 ENVIRONMENT: Indoor Indoor/Outdoor Outdoor/Indoor
 ACTIVITY LEVEL: High Medium Low
 ADJUSTABILITY: Folding Back-height Back angle
 Seat angle COG Q/R wheel
 CONFIGURABILITY: TAB Full arm E/leg rest
 H/Adj arm 8 x 2" S/depth

| NAME | SATURN PORTER WHEELCHAIR | OMEGA ECONOMY |
|--------------|--|--|
| PRODUCT CODE | POR020 | OME018 |
| SIZE | 20" | 18"(NO OPTIONS AVAILABLE) |
| COLOURS | BLACK | BLUE |
| FRAME | NON FOLDING MILD STEEL POWDER COATED REINFORCING PUSH BAR AND BAR BETWEEN FOOTREST HANGERS | POWDER COATED MILD STEEL |
| ARMS | FIXED FULL ARM | REMOVABLE ARMS |
| FOOTREST | FIXED FOOTREST WITH REINFORCING BAR | REMOVABLE |
| CASTORS | 8X2" PU CASTOR ON 2MM STEEL FORK | 8X1" SOLID CASTOR |
| RIMS | 24X1 ^{3/8} MAG RIM WITH 24X1 ^{3/8} PU TYRE - FIXED AXLE | 24X1 ^{3/8} MAG RIM WITH 24X1" SOLID RUBBER TYRE |
| UPHOLSTERY | BLACK REINFORCED NYLON WITH "Z" FOLD POCKET FOR FILES | BLACK REINFORCED NYLON |
| STRENGTH | CHROME PUSHRIMS, STEEL SIDE PANELS, RIGID REINFORCED BOX SHAPE DESIGN | |
| WEIGHT | LOAD CAPACITY - 125KG WEIGHT 24KG | LOAD CAPACITY 110 KG GROSS WEIGHT 20KG |
| PACKING | 110 X 80 X 97 (CM) | 100 X 34 X 85 (CM) |
| ACCESSORIES | STAINLESS STEEL DRIP STAND "T" ROD | |

SATURN ALU-LITE



OMEGA GRANNY (OMEGRA)



| | | | |
|------------------|--|--|--|
| USE: | <input type="checkbox"/> Active | <input checked="" type="checkbox"/> Med-Active | <input type="checkbox"/> Passive |
| ENVIRONMENT: | <input type="checkbox"/> Indoor | <input checked="" type="checkbox"/> Indoor/Outdoor | <input type="checkbox"/> Outdoor/Indoor |
| ACTIVITY LEVEL: | <input type="checkbox"/> High | <input checked="" type="checkbox"/> Medium | <input type="checkbox"/> Low |
| ADJUSTABILITY: | <input checked="" type="checkbox"/> Folding | <input checked="" type="checkbox"/> Back-height | <input type="checkbox"/> Back angle |
| | <input checked="" type="checkbox"/> Seat angle | <input checked="" type="checkbox"/> COG | <input checked="" type="checkbox"/> Q/R wheel |
| CONFIGURABILITY: | <input checked="" type="checkbox"/> TAB | <input type="checkbox"/> Full arm | <input checked="" type="checkbox"/> E/leg rest |
| | <input checked="" type="checkbox"/> H/Adj arm | <input checked="" type="checkbox"/> 8 x 2" | <input checked="" type="checkbox"/> S/depth |

| | | | |
|------------------|---|---|---|
| USE: | <input type="checkbox"/> Active | <input type="checkbox"/> Med-Active | <input checked="" type="checkbox"/> Passive |
| ENVIRONMENT: | <input checked="" type="checkbox"/> Indoor | <input type="checkbox"/> Indoor/Outdoor | <input type="checkbox"/> Outdoor/Indoor |
| ACTIVITY LEVEL: | <input type="checkbox"/> High | <input type="checkbox"/> Medium | <input checked="" type="checkbox"/> Low |
| ADJUSTABILITY: | <input checked="" type="checkbox"/> Folding | <input type="checkbox"/> Back-height | <input type="checkbox"/> Back angle |
| | <input type="checkbox"/> Seat angle | <input type="checkbox"/> COG | <input checked="" type="checkbox"/> Q/R wheel |
| CONFIGURABILITY: | <input type="checkbox"/> TAB | <input type="checkbox"/> Full arm | <input type="checkbox"/> E/leg rest |
| | <input type="checkbox"/> H/Adj arm | <input type="checkbox"/> 8 x 2" | <input type="checkbox"/> S/depth |

| NAME | SATURN ALU-LITE | OMEGA GRANNY |
|--------------|--|---|
| PRODUCT CODE | LW25-77 | OMEGRA |
| SIZE | 14", 16", 18" AND 20" | 16", 18" (NO OPTIONS) |
| COLOURS | METALLIC BLUE, WET BLACK | |
| FRAME | FOLDING ALUMINIUM POWDER COATED | POWDER COATED ALUMINIUM |
| ARMS | FIXED HEIGHT FLIP BACK DESK ARM | FLIP UP DESK ARMS |
| FOOTREST | REMOVABLE | REMOVABLE |
| CASTORS | 8X1" SOLID CASTOR ON ALUMINIUM FORK | 8X1" SOLID CASTOR |
| RIMS | 24X1 ^{3/8} MAG RIM WITH PU TYRE - 12MM QR AXLE | 24X1 ^{3/8} " ALUMINIUM SPOKE RIM WITH 24X1 ^{3/8} " GREY PU TYRE |
| UPHOLSTERY | BLACK REINFORCED NYLON | BLACK REINFORCED NYLON |
| ADJUSTMENT | CENTRE OF GRAVITY, LOW SLUNG, HEIGHT ADJUSTABLE BACKREST WITH PUSH HANDLES | |
| OPTIONS | ELEVATING LEG REST, HEIGHT ADJUSTABLE ARMS, 8X2" CASTOR ON STEEL FORK, ALUMINIUM SPOKE RIM, TENSION ADJUSTABLE BACK UPHOLSTERY, AMPUTEE BRACKET, ANTI-TIP WHEELS | |
| WEIGHT | LOAD CAPACITY 115KG GROSS WEIGHT 14 KG | LOAD CAPACITY - 90KG & GROSS WEIGHT - 20KG |
| PACKING | 100 X 34 X 85 (CM) | 84 X 78 X 32 (CM) |

INVACARE ALU-LITE



PRODUCT CODE: 1517738

Folding Backrest
Reduces overall size when folded for transportation or storage.

Swing away footrest
Footplates fold up as shown and footrests swing away, allowing easier access in or out of the chair.

Foldable cross frame
Lightweight frame folds flat for easier portability.

| | | | |
|------------------|---|---|---|
| USE: | <input type="checkbox"/> Active | <input type="checkbox"/> Med-Active | <input checked="" type="checkbox"/> Passive |
| ENVIRONMENT: | <input checked="" type="checkbox"/> Indoor | <input type="checkbox"/> Indoor/Outdoor | <input type="checkbox"/> Outdoor/Indoor |
| ACTIVITY LEVEL: | <input type="checkbox"/> High | <input type="checkbox"/> Medium | <input checked="" type="checkbox"/> Low |
| ADJUSTABILITY: | <input checked="" type="checkbox"/> Folding | <input type="checkbox"/> Back-height | <input type="checkbox"/> Back angle |
| | <input type="checkbox"/> Seat angle | <input type="checkbox"/> COG | <input type="checkbox"/> Q/R wheel |
| CONFIGURABILITY: | <input type="checkbox"/> TAB | <input type="checkbox"/> Full arm | <input type="checkbox"/> E/leg rest |
| | <input type="checkbox"/> H/Adj arm | <input type="checkbox"/> 8 x 2" | <input type="checkbox"/> S/depth |

Technical data

| | | | | | | | | |
|-------------------|--------------|---------------------|--------|---------|---------|--------|-----|--------------|
| | | | | | | | | |
| Invacare Alu Lite | 405 / 455 mm | 430 mm | 0° | 490 mm | 0° | 450 mm | 80° | 340 - 460 mm |
| | | | | | | | | |
| Invacare Alu Lite | 230 mm | Seat width + 110 mm | 250 mm | 1030 mm | 760 mm | 945 mm | | |
| | | | | | | | | |
| Invacare Alu Lite | 12,5° | 890 mm | 100 kg | 12 kg | 10,2 kg | | | |

BASIK WHEELCHAIR



BASIK ALUMINIUM WHEELCHAIR



USE: Active Med-Active Passive
 ENVIRONMENT: Indoor Indoor/Outdoor Outdoor/Indoor
 ACTIVITY LEVEL: High Medium Low
 ADJUSTABILITY: Folding Back-height Back angle
 Seat angle COG Q/R wheel
 CONFIGURABILITY: TAB Full arm E/leg rest
 H/Adj arm 8 x 2" S/depth

USE: Active Med-Active Passive
 ENVIRONMENT: Indoor Indoor/Outdoor Outdoor/Indoor
 ACTIVITY LEVEL: High Medium Low
 ADJUSTABILITY: Folding Back-height Back angle
 Seat angle COG Q/R wheel
 CONFIGURABILITY: TAB Full arm E/leg rest
 H/Adj arm 8 x 2" S/depth

| NAME | BASIK WHEELCHAIR | BASIK ALUMINIUM WHEELCHAIR |
|--------------|--|---|
| PRODUCT CODE | Yk9031 | ALK903LQB |
| SIZE | 18" | 18" |
| COLOURS | BLACK | |
| FRAME | MILD STEEL POWDER COATED | POWDER COATED ALUMINIUM |
| ARMS | FLIP BACK DESK ARMS | FLIP UP DESK ARMS |
| FOOTREST | REMOVABLE WITH PLASTIC FOOTPLATE | REMOVABLE |
| CASTORS | 8X1" SOLID CASTOR | 8X1" SOLID CASTOR |
| RIMS | 24X1" SPOKE RIM WITH PLASTIC PUSHRIM 24X 1" SOLID RUBBER TYRE | 24X1 ^{3/8} " ALUMINIUM SPOKE RIM WITH 24X 1 ^{3/8} " GREY PNEUMATIC |
| UPHOLSTERY | BLACK NYLON | BLACK NYLON |
| WEIGHT | LOAD CAPACITY 100KG GROSS WEIGHT 19KG | LOAD CAPACITY 85KG GROSS WEIGHT 13KG |
| PACKING | 100 X 34 X 85 (CM) | 84 X 78 X 32 (CM) |

BASIK TRANSIT



BASIK HYBRID GRANNY



- Wheelchair/Transport chair two way use
- Aluminium chair frame
- Flip back desk arm
- Detachable mag wheel with PU tyre
- 8" rear wheel & 8" front castor
- Tool-free adjustable swing away footrest
- Aluminium front fork with 2 adjustable positions

| | | | |
|------------------|---|---|---|
| USE: | <input type="checkbox"/> Active | <input type="checkbox"/> Med-Active | <input checked="" type="checkbox"/> Passive |
| ENVIRONMENT: | <input checked="" type="checkbox"/> Indoor | <input type="checkbox"/> Indoor/Outdoor | <input type="checkbox"/> Outdoor/Indoor |
| ACTIVITY LEVEL: | <input type="checkbox"/> High | <input type="checkbox"/> Medium | <input checked="" type="checkbox"/> Low |
| ADJUSTABILITY: | <input checked="" type="checkbox"/> Folding | <input type="checkbox"/> Back-height | <input type="checkbox"/> Back angle |
| | <input type="checkbox"/> Seat angle | <input type="checkbox"/> COG | <input type="checkbox"/> Q/R wheel |
| CONFIGURABILITY: | <input type="checkbox"/> TAB | <input type="checkbox"/> Full arm | <input type="checkbox"/> E/leg rest |
| | <input type="checkbox"/> H/Adj arm | <input type="checkbox"/> 8 x 2" | <input type="checkbox"/> S/depth |

| | | | |
|------------------|---|---|---|
| USE: | <input type="checkbox"/> Active | <input type="checkbox"/> Med-Active | <input checked="" type="checkbox"/> Passive |
| ENVIRONMENT: | <input checked="" type="checkbox"/> Indoor | <input type="checkbox"/> Indoor/Outdoor | <input type="checkbox"/> Outdoor/Indoor |
| ACTIVITY LEVEL: | <input type="checkbox"/> High | <input type="checkbox"/> Medium | <input checked="" type="checkbox"/> Low |
| ADJUSTABILITY: | <input checked="" type="checkbox"/> Folding | <input type="checkbox"/> Back-height | <input type="checkbox"/> Back angle |
| | <input type="checkbox"/> Seat angle | <input type="checkbox"/> COG | <input checked="" type="checkbox"/> Q/R wheel |
| CONFIGURABILITY: | <input type="checkbox"/> TAB | <input type="checkbox"/> Full arm | <input type="checkbox"/> E/leg rest |
| | <input type="checkbox"/> H/Adj arm | <input type="checkbox"/> 8 x 2" | <input type="checkbox"/> S/depth |

| NAME | BASIK TRANSIT WHEELCHAIR | BASIK HYBRID GRANNY |
|--------------|--|--|
| PRODUCT CODE | YK9092/3 | Yk9061 |
| SIZE | 17" & 19" | 18" |
| COLOURS | BLACK | SILVER |
| FRAME | MILD STEEL & ALUMINIUM POWDER COATED | ALUMINIUM CHAIR FRAME |
| ARMS | FIXED FULL ARMS | FLIP BACK DESK ARM |
| FOOTREST | REMOVABLE | TOOL-FREE ADJUSTABLE SWING AWAY FOOTREST |
| CASTORS | FOUR 8X1" SOLID CASTOR | 24" REAR WHEEL & 8" FRONT CASTOR |
| RIMS | | FIXED MAG WHEEL WITH PU TYRE |
| UPHOLSTERY | BLACK NYLON | BLACK NYLON |
| WEIGHT | LOAD CAPACITY – 90KG GROSS WEIGHT – 11KG | LOAD CAPACITY 90KG & WEIGHT 16KG |
| PACKING | 76 X 26 X 60 (CM) | |
| USE | LIGHT INDOOR / SMOOTH SURFACE USE | LIGHT INDOOR / SMOOTH SURFACE USE |