

ANGLED GRAB RAILS



PRODUCT CODE: DOGLEG

PRODUCT CODE: 45GRAB

HINGED GRAB RAIL



*Illustration of folding



PRODUCT CODE: QSE-088Q

GRAB RAILS CHROME OR STAINLESS STEEL



PRODUCT CODE: APC211CHROME



PRODUCT CODE: SG088S/STEEL



STILEO SUCTION GRAB RAIL



NAME	DOG-LEG S/STEEL & 45° CHROME GRAB RAIL	HINGE OR FOLD-UP GRAB RAIL	SUCTION UNIT – HOSPITAL USE	SUCTION GRAB RAIL –INVACARE H171 STILEO
PRODUCT CODE	DOGLEG	QSE-088Q	APC211CHROME	H171
SIZE HEIGHT	EACH SECTION MEASURES 3CM 70CM – DOG-LEG	750 MM	12" (30CM), 18" (46CM), 24" (61CM) OPTION OF 35MM OR	40CM
MATERIAL	35CM - 45° DOG-LEG – STAINLESS STEEL 45° CHROMED MILD STEEL WITH GNARLED GRIP	STAINLESS STEEL	70MM IN CHROME 35MM ONLY IN S/STEEL CHROMED IN MILD STEEL WITH GNARLED GRIP S/STEEL ROUGH POLISH	MOULDED PLASTIC WITH GRIP KNURLING
SUPPLIED	SCREWS (PLUG NOT INCLUDED)	SCREWS (PLUG NOT INCLUDED)	SCREWS (PLUG NOT INCLUDED)	
COLOUR				BLUE – WHITE
LOAD CAPACITY	100KG	100KG	100KG	100KG
WEIGHT		WEIGHT – 3.2KG EACH GROSS WEIGHT – 4.5KG		WEIGHT – 0.5KG EACH
FUNCTION				VACUUM SUCTION WITH GREEN / RED INDICATOR LEVELS
PACKING		800X320X120 (MM)		

TOILET SEAT RAISERS



Invacare H304

- Toilet seat raiser made of soft, skin friendly material
- Raises the toilet set by 9cm
- Fits on the toilet easily and safely, no assembly required
- Large hygiene recess at the front for easy personal hygiene

Invacare H301

- Toilet seat raiser
- Very easy to attach to the toilet
- Smooth surface for easy cleaning
- Large hygiene recess (24 x 29 cm)



Aquatec® 90 - AT90

- Toilet seat raiser
- Raises the toilet seat by 10cm
- Attached safely to the toilet with two side clamps
- Large turning knobs for quick attachment and detachment
- Particularly high load capacity: 226kg

TOILET SEAT RAISERS



PRODUCT CODE: AT900



PRODUCT CODE: 1302RTS

Aquatec® 900

- Securely attached toilet seat raiser
- Individual adjustment to the user at any time, no tools required:
 - Three different seat heights
 - Three different angles, lower at the front
- Fold-way, ergonomic armrests for additional support when sitting down and standing up
- Available with or without armrests
- Seat can be folded up in all setting positions – ease of cleaning

Invacare® Raised Toilet Seat Model no. 1301RTS without arms

Model no. 1302RTS with arms (shown)

- Comfortable, contoured seat is made of easy-to-clean, durable moulded composite
- Fits on most round bowl configurations
- Clamp-on locking mechanism tightens easily with one hand
- New adjustable lock arms (13302RTS)

Technical data

Order Number	Invacare H304	Invacare H301	Aquatec® 90
Toilet seat raiser	1470982		
Toilet seat raiser with lid		1522262	8.01.001
Toilet seat raiser without lid			8.00.001

Dimensions			
Seat width	38.5cm (A)	36.5 cm (A)	36.5 cm (A)
Seat depth	43cm (B)	38.5 cm (B)	39.5 cm (B)
Seat height	9 cm (C)	-	-
Seat height front	-	12.5 (cm) (C)	10.5 cm (C)
Seat height rear	-	13.5 cm (D)	11.5 cm (D)
Total width	40cm	38cm	-
Total width (with lid)	-	-	38 cm (E)
Total width (without lid)	-	-	37.5 cm (F)
Total depth	43.5cm (F)	40 cm (F)	41 cm (F)
Total height	13.5 cm (G)	18.5 cm (G)	-
Total height (with lid)	-	-	13.5 cm
Total height (without lid)	-	-	12.5 cm (G)
Width of opening	22.5 cm (H)	19 cm (H)	21 cm (H)
Depth of opening	27.5 cm (I)	25 cm (I)	27 cm (I)

More Information			
Weight	2.2kg	0.9kg	-
Weight (with lid)	-	-	1.6kg
Weight (without lid)	-	-	1.1kg
Load capacity	125kg	120kg	225kg
Colour	Grey	White	White