

STATIO BED CRADLE



ARESTE BACK SUPPORT



<b>NAME</b>	INVACARE STATIO BED CRADLE	INVACARE ARESTE BACKREST LIFT ASSIST
<b>PRODUCT CODE</b>	L871	L110
<b>FRAME</b>	POWDER COATED MILD STEEL	POWDER COATED MILD STEEL WITH BREATHABLE MESH
<b>BASE</b>	54L X 32 W X 34 H (CM)	510X790 (MM)
<b>CARRY CAPACITY</b>		5 ADJUSTABLE POSITIONS WEIGHT – 100KG GROSS WEIGHT – 4.5KG
<b>PACKING</b>		600X660X125 (MM)

PONA FOLDING TABLE



ADEO TOILET FRAME



PRODUCT CODE: 970790

PRODUCT CODE: C407A

SUCTION CLAW & REACHER



UNDER MATTRESS COTSIDE



\* Chrome Finish

<b>NAME</b>	REACHER WITH SUCTION	<b>NAME</b>	UNDER MATTRESS COT SIDES
<b>PRODUCT CODE</b>	APC40130	<b>PRODUCT CODE</b>	COT001
<b>REACH</b>	REACHER WITH CLAW (10209)	<b>FRAME</b>	POLISHED CHROME MILD STEEL
<b>PACKING</b>	ADJUSTABLE SLATTED MOULDED PLASTIC 1 PER BOX / PACKET	<b>ADJUSTABLE</b>	HEIGHT – 220 TO 550 (MM) LENGTH – 1,180 TO 1,800 (MM)
<b>REACH</b>	80CM	<b>TO FIT</b>	SINGLE BED
<b>LOAD CAPACITY</b>	2KG	<b>WEIGHT</b>	GROSS WEIGHT – 2.5KG
		<b>PACKING</b>	1,210X560X60 (MM)

TRACTION UNITS

OVERDOOR



PRODUCT CODE: 24300

MATTRESS  
(can also be used for neck)



PRODUCT CODE: 24330

OVERBED TABLE/TILT



PRODUCT CODE: APC10210



PRODUCT CODE: APC10213

MOBILE IV STAND



<b>NAME</b>	OVERBED TABLE (WITH TILT)
<b>PRODUCT CODE</b>	APC10210 & APC10213
<b>FRAME</b>	MILD STEEL POWDER COATED
<b>CASTOR</b>	2" PLASTIC
<b>ADJUSTABILITY</b>	60-114CM 76-106CM (TILT)
<b>MECHANISM</b>	SPRING LOADED MANUAL (TILT)
<b>SURFACE</b>	DARK WOOD MELAMINE LAMINATE 38X80CM 38-60CM & 38-20 CM TILT) – SPLIT SURFACE
<b>COLOUR</b>	SILVER VEIN
<b>LOAD CAPACITY</b>	6KG ON SURFACE
<b>GROSS WEIGHT</b>	11.2KG 10.2 KG (TILT)

<b>NAME</b>	MOBILE DRIP STAND
<b>PRODUCT CODE</b>	IV
<b>BASE</b>	MOULDED PLASTIC WITH PLASTIC 2" CASTOR EACH LEG EXTENDS – 25CM
<b>T ROD</b>	CHROME POLISH MILD STEEL
<b>WEIGHT</b>	ADJUSTABILITY – 1,050 TO 1,050 (MM) GROSS WEIGHT – 4KG EACH
<b>PACKING</b>	4 PER BOX – 1,490X140X570(MM)

MONKEY CHAIN FREESTANDING



INVACARE PORTO LIFT OUT CHAIR



**Relax position**  
Different relax positions will allow the same comfort in any recliner angle chosen.



**Normal position**  
Different fabrics will allow a perfect fit in all environments.

PRODUCT CODE: PORTO

The Porto chair provides a high level of comfort thanks to the continuing of contours between the leg rest, the seat and the lumbar support which maximizes good positioning. Porto has a contemporary design and can be easily moved from one place to another, thanks to the small back castors which are simply dismantled by taking out a screw.

- One motor
- Easy-to-use remote control
- Side pocket
- Easy to use and mount

<b>NAME</b>	MONKEY CHAIN – FREE STANDING
<b>PRODUCT CODE</b>	SMA11000
<b>FRAME</b>	POWDER COATED MILD STEEL
<b>BASE</b>	102LX76WX97H (CM) ON ADJUSTABLE FEET FOR STABILITY
<b>TOP BAR</b>	HEIGHT 163 (CM) REACH 91 (CM) MOUNTING CLAMP FOR 2 TO 7.6 CM) TUBING
<b>CARRY CAPACITY</b>	WEIGHT – 110KG GROSS WEIGHT – 21KG
<b>PACKING</b>	1060X770X120 (MM)

**Positions**

**Colour**

**INVACARE BEDPAN AND URINAL BOTTLE**



**BASIK BEDPAN AND URINAL BOTTLE**



ALK001  
Urinal

PRODUCT CODE: ALK001



ALK002  
Urinal

PRODUCT CODE: ALK002

<b>NAME</b>	INVACARE BEDPAN AND URINAL BOTTLE
<b>PRODUCT CODE</b>	BEDPAN 970765   URINAL BOTTLE 970742
<b>FRAME</b>	MOULDED PLASTIC
<b>DIMENSIONS</b>	BEDPAN: TOTAL TO END OF HANDLE – 470MM 380MM TO END OF OPENING OPENING OF 200MM URINAL BOTTLE: TOTAL OF 300MM NECK – 100MM



ALK003  
Urinal

PRODUCT CODE: ALK003