

# INVACARE Action Family: The changing face of Mobility



**The Action Family product offer consists of:**  
 Action2 - Standard aluminium wheelchair  
 Action3 Junior - Paediatric wheelchair from age 3 and up

Action3 - Lightweight wheelchair, single crossbar, highly configurable with a huge range of options, including:

- Hemiplegic version
- One arm drive version
- Double hand rim version
- Comfort version
- and more

Action4 - Lightweight wheelchair, double crossbar, highly configurable with huge range of options, including:

- Hemiplegic version
  - One arm drive version
  - Double hand rim version
  - Comfort version
- Heavy-duty option up to 160kg user weight

Cross Compatibility table:	Action <sup>2</sup>	Action <sup>3</sup>
Leg rest	x	x
Footplates	x	x
Armrest	x	x
Arm pads	x	x
Side guards	x	x
Seat upholstery	x	x
Backrest canes	x	x
Backrest upholstery		x
Push handles	x	x
Brakes	x	x
Fork		x
Castor wheels	x	x
Frame	x	x
User weight	125KG	125KG

Action<sup>3</sup> One arm drive version



Action<sup>3</sup> Comfort version



Action<sup>3</sup> Hemiplegic version



INVACARE ACTION 1



**Legrest**  
Different types of legrest are available: amputee support , elevating legrest 0° to 80° and fixed legrest , to meet every user needs



**Headrest**  
Can be adjusted in both height and depth to the backbrace to bring more comfort to the user.



**Front castor**  
2 types of front castor are available: 8 x 1.3/8" solid and 8" x 2 solid

INVACARE ACTION 2



Features



**Swing in and out footrest**  
The 80 degree footrest swings in and out for easy transfers in or out of the chair



**Adjustable center of gravity**  
The newly-designed rear wheel pane can be adjusted by 5 cm.

PRODUCT CODE: ACT1

USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input checked="" type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input checked="" type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input checked="" type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input checked="" type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

PRODUCT CODE: ACT2

USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input checked="" type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input type="checkbox"/> Medium	<input checked="" type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input checked="" type="checkbox"/> COG	<input checked="" type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

Self propulsion: 300 / 300 / 400 / 400 / 400 mm	420 mm	31"	400 / 500 mm	31"	310 mm	80° / 0° to 80°
Self propulsion: 300 / 400 mm	320 mm	Seat width = 210 / 300 mm**	210 / 300 mm**	3048 / 1068 mm**	310 / 320 mm**	800 / 910 mm
Self propulsion: 127 kg	16.5 / 18.8 / 16.0 kg**	16.5 / 11.8 / 12.9 kg**				

Technical data							
36-40.5-41-40.5-48 cm	40.5 / 45 cm	46.5-51 cm	61 cm	20-36 cm	35- 60 cm	80° / 0- 80°	9° / 0° - 30°
100+ 200 cm	112 cm	82 / 88.2 cm	41 cm	15 kg	max 120 kg	Crab wheel	

## INVACARE ACTION 3



## INVACARE ACTION 4



PRODUCT CODE: ACT3

USE:	<input checked="" type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input checked="" type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input checked="" type="checkbox"/> COG	<input checked="" type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input checked="" type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input checked="" type="checkbox"/> H/Adj arm	<input checked="" type="checkbox"/> 8 x 2"	<input checked="" type="checkbox"/> S/depth

PRODUCT CODE: ACT4/HD

USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input checked="" type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input checked="" type="checkbox"/> COG	<input checked="" type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input checked="" type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input checked="" type="checkbox"/> H/Adj arm	<input checked="" type="checkbox"/> 8 x 2"	<input checked="" type="checkbox"/> S/depth

38-46.5, 41, 43.5, 46, 48.5, 46, 50.5 cm	41, 46 cm	38.5, 41, 43.5, 46, 48.5, 51 cm	41, 45.5, 51 cm	30-33-36 cm, 24 cm, 19-23.5 cm	31-48 cm, 34-47 cm	87° (7°-87°)
17, 17.5 (17°-23°), 18°-23°	Seat width = 23 cm	81.5-102 cm	100 cm	72-85 cm	38.5 cm	14.2 kg
			* Our wheelchair comply above ISO standards and with ISO 13242, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.			
144 (12) kg	Transport weight = 74 kg	Castor locked (front, armrest, footrest and rear wheel)	ISO 7176-19*			

38, 40.5, 43, 45.5, 48, 50.5, 53.5, 46.5 cm	41, 46 and 50 cm	38.5, 41, 43.5, 46, 48.5, 51 cm	38-48 cm, 38-51 cm**	20-23-26 cm, 34 cm, 19-23.5 cm	31-48 cm	0°-4°
7°-10°	Seat width = 21 cm	80-103 cm	105-111 cm	81-90 cm	50.5 cm	13.2 kg
			* January 2018 ** Reclining backrest *** Our wheelchair comply above ISO standards and with ISO 13242, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.			
132 (10) kg	Transport Weight 8.3 kg	Castor locked ISO 7176-19***	ISO 7176-19***			

INVACARE ACTION 5



The Invacare Action5 has been designed for the more active user looking for an affordable chair with good drive performance and the ability to fold easily for handling, storage and transportation purposes. With its two patents and lightweight construction, the new Action5 gives the user a chair they can rely on. Key features include:

- Unique folding mechanism
- Stability of a rigid chair
- Quick step-less adjustability
- Highly configurable

USE:	<input checked="" type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input checked="" type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input checked="" type="checkbox"/> COG	<input checked="" type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input checked="" type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input checked="" type="checkbox"/> H/Adj arm	<input checked="" type="checkbox"/> 8 x 2"	<input checked="" type="checkbox"/> S/depth

PRODUCT CODE: SJC0110

TECHNICAL DATA

<i>Invacare Action<sup>5</sup></i>	305 / 330 / 355 / 380 / 405 / 430 / 455 / 480 / 505 / 530 mm	350 / 375 / 400 / 425 / 450 / 475 / 500 mm	0 - 9°	Front 360-535 mm in steps of 25 mm Rear: 360-510 mm in steps 25 mm	Fixed 90° Or stepless 90° -15°/+15°	Fixed 300 / 325 / 350 / 400 / 425 / 450 / 475 / 500 mm Adjustable 360-400 and 410-490 mm	Fixed front 80° Swing away legrest 70° / 80° & 90° Elevating legrest 0-80°
<i>Invacare Action<sup>5</sup></i>	Swing away frame: 200-510 mm Fixed frame: 300-530 mm	186-227 mm 186-266 mm 205-280 mm	SW +160 mm To SW+ 250 mm	320 -380 mm	845 -1205 mm	670 -860 mm	701 -1108 mm
<i>Invacare Action<sup>5</sup></i>	7°	1050 -1380 mm	100 & 130 kg	From 12,3 kg	From 8,3 kg		