



Yes, you can.®

küschall® K-SERIES



Manual Wheelchair aluminium, active, rigid frame, folding backrest and angle adjustable

K4 frame : hangers 75° aluminium, carbon or titane
 Air Lite frame : hangers 90°, aluminium or titane frame
 New adducted frame for K4 and Air Lite aluminium
 Weight : from 7,7 kg* full (7,5 kg carbon frame)
 Max user weight user : 120 kg

Date :
 User name:
 Dealer stamp:

* Configuration 7,7kg:SB40cm, ST40 cm,RH 37,5 cm, backrest tube titane, footrest titane, wheel light fibre, wheel 3° supersport, camber tube carbon

(nc 9999) : not compatible with option 9999

Euros
Public Price
HT
2 417,06

DDA0067	küschall® K-Series	2 417,06
---------	---------------------------	----------

Seat width, seat heights, backrest height and UL length could vary with +/-10 mm.

Please make your selection with a cross in the box : Standard feature Free of Charge

Select your Pac-K Series and beneficiate a price advantage :

⚠ Pac-K Series (The weights indicated must be added or subtract from 9,9 kg, weight for stanlyard configuration in 400x400 mm, wheel Starec)/ if you select a Pack-K, please do not tick the option individually - you can choose only one pack

DDA2000	Pac-K Series Carbon : Carbon frame 75° - Carbon footplate - Carbon axle +150g	<input type="checkbox"/>	796,21
DDA2001	Pac-K Series Titane : Titanium frame - Titanium backrest - Titanium footrest -225g	<input type="checkbox"/>	843,60
	Titanium frame 75° <input type="checkbox"/>		Titanium frame 90° <input type="radio"/>
DDA2002	Pac-K Series 1 : Titanium backrest - Titanium footrest - Carbon axle -305g	<input type="checkbox"/>	426,54
DDA2003	Pac-K Series 2 : Titanium backrest - Titanium footrest - Carbon axle - Seat upholstery light - Backrest upholstery light -720g	<input type="checkbox"/>	559,24
DDA2004	Pac-K Series 3 : Titanium backrest - Titanium footrest - Carbon axle - Seat upholstery light - Backrest upholstery light -Spinergy LXX rear wheels - Speedrun tyre - Active brakes	<input type="checkbox"/>	938,39
DDA2008	Brave pack : Bottle holder - Mobile phone holder - Mountain bike rear wheels 24" -1470g	<input type="checkbox"/>	426,54

SEAT TYPE

FRAME (for tapered frame: + 10 mm of frame lenght)

DDA0284	K4, 75° to the ground Aluminium (nc ST 475 AND 525 MM°)	<input type="checkbox"/>	<input checked="" type="checkbox"/>
DDA0285	K4, 75° to the ground Aluminium adducted (nc ST 475 et 525 mm)	<input type="checkbox"/>	<input type="radio"/>
DDA0286	K4, 75° to the ground Aluminium +50 mm (nc ST 375, 400 mm)	<input type="checkbox"/>	102,37
DDA0287	K4, 75° to the ground Titanium (nc ST 375, 475, 500, 525 mm, nc SHv 480 mm, nc upmounted footrest nc one piece swing away footrest, nc frame colours, nc suspensions, nc rear wheel extension)	<input type="checkbox"/>	654,03
DDA0288	K4, 75° to the ground Carbon (nc ST 475 et 500 mm / nc UL < 340 mm / nc SHv 480 mm, nc Airlite Frame, hangers 90°, nc rear wheel extension) -200g	<input type="checkbox"/>	748,82
DDA0289	AirLite, 90° to the ground Aluminium (nc wheel 5") (nc ST 475 et 525 mm)	<input type="checkbox"/>	156,40
DDA0290	AirLite, 90° to the ground Aluminium adducted (nc wheel 5")(nc ST 475 et 525 mm)	<input type="checkbox"/>	156,40



Adducted frame




DDA0291	AirLite, 90° to the ground Aluminium +50 mm (nc wheel 5") (nc ST 475 et 525 mm)	<input type="checkbox"/>	156,40
DDA0292	AirLite, 90° to the ground Titanium (nc wheel 5")(nc ST 475 et 525 mm)	<input type="checkbox"/>	654,03




SEAT WIDTH(SB) (SB 340, 440, 460 mm)

DDA0124	340 mm	<input type="checkbox"/>	DDA0130	400 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0126	360 mm	<input type="checkbox"/>	DDA0132	420 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0128	380 mm	<input type="checkbox"/>	DDA0134	440 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0136	460 mm	<input type="checkbox"/>			<input type="checkbox"/>	<input type="radio"/>
DDA0138	480 mm (nc Alcantara upholstery)	<input type="checkbox"/>			<input type="checkbox"/>	<input type="radio"/>




SEAT DEPTH (ST) (475 to 525 mm only with long frame)

DDA0414	375 mm	<input type="checkbox"/>	DDA0423	450 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0417	400 mm	<input type="checkbox"/>	DDA0426	475 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0420	425 mm	<input type="checkbox"/>	DDA0429	500 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0432	525 mm	<input type="checkbox"/>			<input type="checkbox"/>	<input type="radio"/>



SEAT POSITION, possibilities of adjustments

DDA0461	140 mm Position 1	<input type="checkbox"/>	DDA0462	125 mm Position 2	<input type="checkbox"/>	<input type="radio"/>
DDA0463				110 mm Position 3	<input type="checkbox"/>	<input checked="" type="radio"/>
DDA0464	95 mm Position 4	<input type="checkbox"/>	DDA0465	80 mm Position 5	<input type="checkbox"/>	<input type="radio"/>



BACK HEIGHT(RH)

Short


DDA0514	375 mm	<input type="checkbox"/>			<input type="checkbox"/>	<input checked="" type="radio"/>
DDA0548	270 mm	<input type="checkbox"/>	DDA0506	330 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0580	285 mm	<input type="checkbox"/>	DDA0509	345 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0550	300 mm	<input type="checkbox"/>	DDA0512	360 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0503	315 mm	<input type="checkbox"/>			<input type="checkbox"/>	<input type="radio"/>

Long

DDA0516	390 mm	<input type="checkbox"/>			<input type="checkbox"/>	<input type="radio"/>
DDA0518	405 mm	<input type="checkbox"/>	DDA0524	450 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0520	420 mm	<input type="checkbox"/>	DDA0526	465 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0522	435 mm	<input type="checkbox"/>	DDA0528	480 mm	<input type="checkbox"/>	<input type="radio"/>

FRONT SEAT-TO-FLOOR HEIGHT (SHv) (Refer to the appendix page 31)

DDA0847	450 mm only with standard frame	<input type="checkbox"/>	<input type="radio"/>
DDA0848	460 mm with standard frame	<input type="checkbox"/>	<input type="radio"/>
DDA0849	470 mm with standard frame	<input type="checkbox"/>	<input type="radio"/>
DDA0850	480 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0851	490 mm	<input type="checkbox"/>	<input checked="" type="radio"/>
DDA0852	500 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0853	510 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0854	520 mm only with long frame (adjustments until 540 mm with wheel 26" in case of special built)	<input type="checkbox"/>	<input type="radio"/>



REAR SEAT-TO-FLOOR HEIGHT(SHh) (with suspensions, SHh delivered + 30 mm / requested height) (SHh 470, 480 et 490 nc suspensions)

DDA0948	430 mm	<input type="checkbox"/>	<input checked="" type="radio"/>			
DDA0943	380 mm (nc 25", 26")	<input type="checkbox"/>	DDA0950	450 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0944	390 mm (nc 26")	<input type="checkbox"/>	DDA0951	460 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0945	400 mm	<input type="checkbox"/>	DDA0952	470 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0946	410 mm	<input type="checkbox"/>	DDA0953	480 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0947	420 mm	<input type="checkbox"/>	DDA0954	490 mm	<input type="checkbox"/>	<input type="radio"/>
DDA0949	440 mm	<input type="checkbox"/>			<input type="checkbox"/>	<input type="radio"/>


Type de fixations du module d'assise

SHH	Taille des roues arrière		
	24"	25"	26"
380 mm	S	-	-
390 mm	S	S	-
400 mm	S	S	S
410 mm	S	S	S
420 mm	M	S	S
430 mm	M	M	S
440 mm	M	M	M
450 mm	M	M	M
460 mm	L	M	M
470 mm	L	L	M
480 mm	L	L	L
490 mm	L	L	L

KNEE-TO-HEEL LENGTH
(K4 : UL =< SHv - 30 mm / Airlite : UL =< Shv - 40 mm)

UL from 290 to 330 mm : upmounted footrest (nc protection frame, frame carbon/titane)

DDA1670	290 mm	<input type="checkbox"/>	DDA1673	320 mm	<input type="checkbox"/>	<input type="radio"/>
DDA1671	300 mm	<input type="checkbox"/>	DDA1674	330 mm	<input type="checkbox"/>	<input type="radio"/>
DDA1672	310 mm	<input type="checkbox"/>			<input type="checkbox"/>	<input type="radio"/>



UL from 340 to 480 mm : standard footrest (except for long frame and Airlite : if UL from 340 to 400 mm, upmounted footrest)

DDA1684	430 mm	<input type="checkbox"/>	<input checked="" type="radio"/>			
DDA1675	340 mm	<input type="checkbox"/>	DDA1682	410 mm	<input type="checkbox"/>	<input type="radio"/>
DDA1676	350 mm	<input type="checkbox"/>	DDA1683	420 mm	<input type="checkbox"/>	<input type="radio"/>
DDA1677	360 mm	<input type="checkbox"/>	DDA1685	440 mm	<input type="checkbox"/>	<input type="radio"/>
DDA1678	370 mm	<input type="checkbox"/>	DDA1686	450 mm	<input type="checkbox"/>	<input type="radio"/>
DDA1679	380 mm	<input type="checkbox"/>	DDA1687	460 mm	<input type="checkbox"/>	<input type="radio"/>
DDA1680	390 mm	<input type="checkbox"/>	DDA1688	470 mm	<input type="checkbox"/>	<input type="radio"/>
DDA1681	400 mm	<input type="checkbox"/>	DDA1689	480 mm	<input type="checkbox"/>	<input type="radio"/>

BACKREST

BACKREST MATERIAL

DDA0730	Aluminium	<input type="checkbox"/>	◆
DDA0731	Titanium = - 80 g (with aluminium push handles)	<input type="checkbox"/>	246,45



0731
0722

BACK UPHOLSTERY


DDA0721	Light back upholstery (hand made) (nc push handles height adjustable)	-355g	<input type="checkbox"/>	127,96
DDA0722	Tension velcro adjustable back Black/grey		<input type="checkbox"/>	◆
DDA0735	Aero back, tension adjustable		<input type="checkbox"/>	184,83
DDA0772	Without backrest upholstery (nc: titan backrest, long, court and externals adjustable height push handles)		<input type="checkbox"/>	○

0721



PUSH HANDLES

DDA0750	Without push handles	<input type="checkbox"/>	◆
DDA0751	Push handles fixed, short	<input type="checkbox"/>	80,57
DDA0702	Push handles fixed, long	<input type="checkbox"/>	85,31
DDA0703	Externals Height adjustable, rear set (nc tubular armrests)	<input type="checkbox"/>	151,66
DDA0733	Height adjustable integrated (only for option "without backrest upholstery)	<input type="checkbox"/>	132,70
DDA0704	Foldables	<input type="checkbox"/>	151,66



0703

BACKREST ANGLE (to be measured from the seat)

DDA0707	74°	<input type="checkbox"/>	○
DDA0718	78°	<input type="checkbox"/>	○
DDA0709	82°	<input type="checkbox"/>	◆
DDA0719	86°	<input type="checkbox"/>	○
DDA0711	90°	<input type="checkbox"/>	○



SEAT UPHOLSTERY

DDA0610	Standard	<input type="checkbox"/>	◆	
DDA0615	Toile Air Light	-60g	<input type="checkbox"/>	47,39
DDA0611	Tension adjustable by velcro	<input type="checkbox"/>	66,35	

0615



SEAT UPHOLSTERY BAGS - only with standard upholstery

DDA0616	Small bag	<input type="checkbox"/>	42,65
DDA0617	Big bag	<input type="checkbox"/>	42,65



FOOTREST

FOOTREST TUBE (nc footplate)

DDA1508	Aluminium	<input type="checkbox"/>	◆	
DDA1509	Titane	-145g	<input type="checkbox"/>	175,36
DDA1513	Footrest cover (on standard footrest)	<input type="checkbox"/>	81,52	

1509



FOOTPLATE

DDA1570	Carbon footplate, angle adjustable	<input type="checkbox"/>	189,57
DDA1515	Footplate, flip-up to side (nc adducted frame, carbon/Titane)	<input type="checkbox"/>	327,96

1510



SIDE GUARDS / CLOTHES GUARDS / MUD GUARDS


CARBON CLOTHES GUARD

DDA1401	Adjustable and foldable	<input type="checkbox"/>	◆
---------	-------------------------	--------------------------	---

1401 / 1408

MUD GUARD (CARBON): With suspension, distance between tyre and mud guard is 40 mm

DDA1408	Removable	<input type="checkbox"/>	232,23
---------	-----------	--------------------------	--------



1416

SIDE GUARD (WITH CLOTHES GUARD) (nc SHh < 420 mm)

DDA1416	Siderest, height adjustable, short	<input type="checkbox"/>	199,05
---------	------------------------------------	--------------------------	--------

ARMREST (nc with rear wheel extension in combination with SHh ≤ 410 mm)

DDA1411	Tubular armrest, height adjustable, removable, swing-away, comfort padded	<input type="checkbox"/>	208,53
---------	---	--------------------------	--------

1411

std-262g

SEAT OPTIONS

SEAT CUSHION

DDA0600	Standard H: 80 mm	<input type="checkbox"/>	73,93
---------	-------------------	--------------------------	-------



0606

SEAT BAG

DDA0606	Seat bag Alcantara	<input type="checkbox"/>	59,72
---------	--------------------	--------------------------	-------

FRAME

FRAME COLOUR (nc frame carbon/Titane)

DDA0300	Glossy black	<input type="checkbox"/>	DDA0350	Purple fire	<input type="checkbox"/>	<input type="radio"/>
DDA0358	Smokey grey	<input type="checkbox"/>	DDA0302	Royal blue	<input type="checkbox"/>	<input type="radio"/>
DDA0356	Ivory white	<input type="checkbox"/>	DDA0315	Silver grey	<input type="checkbox"/>	<input type="radio"/>
DDA0303	Royal red	<input type="checkbox"/>	DDA0314	Black matt	<input type="checkbox"/>	<input type="radio"/>
DDA0354	Night shimmer	<input type="checkbox"/>	DDA0366	Dynamite Orange	<input type="checkbox"/>	<input type="radio"/>
DDA0359	Flashy yellow	<input type="checkbox"/>	DDA0367	Atomic Purple	<input type="checkbox"/>	<input type="radio"/>
DDA0365	Mercury Blue	<input type="checkbox"/>	DDA0368	Nitro Pink	<input type="checkbox"/>	<input type="radio"/>
DDA0310	Mirror natural	<input type="checkbox"/>				218,01

included)

DDA0330	techno-anthracite	<input type="checkbox"/>	DDA0334	b-blue	<input type="checkbox"/>	322,27
DDA0336	arctic blue	<input type="checkbox"/>	DDA0329	comet red	<input type="checkbox"/>	322,27
DDA0338	chrome	<input type="checkbox"/>	DDA0337	copper glaze	<input type="checkbox"/>	322,27
DDA0339	Champagne	<input type="checkbox"/>	DDA0340	Olive	<input type="checkbox"/>	322,27

SUSPENSION for SHh 380 -460 mm- user weight (nc seat width > 480 mm)

DDA0630	40 - 65 kg	<input type="checkbox"/>				228,44
DDA0631	65 - 85 kg	<input type="checkbox"/>				228,44
DDA0632	85 - 100 kg	<input type="checkbox"/>				228,44

FRONT WHEEL

DDA0800	3" Supersport	<input type="checkbox"/>				<input type="radio"/>
DDA0801	4" Everyday	<input type="checkbox"/>				<input checked="" type="radio"/>
DDA0860	4" Comfort wheel black rubber	<input type="checkbox"/>				<input type="radio"/>
	4" Starec :					
DDA0803	techno-anthracite	<input type="checkbox"/>	DDA0808	comet red	<input type="checkbox"/>	54,98
DDA0831	artic blue	<input type="checkbox"/>	DDA0809	black	<input type="checkbox"/>	54,98
DDA0861	chrome	<input type="checkbox"/>	DDA0830	copper glaze	<input type="checkbox"/>	54,98
DDA0825	Champagne	<input type="checkbox"/>	DDA0826	Olive	<input type="checkbox"/>	54,98
DDA0807	b-blue	<input type="checkbox"/>				54,98
DDA0824	Skater wheel	<input type="checkbox"/>				19,91
DDA0810	5" Light wheel	<input type="checkbox"/>				<input type="radio"/>
DDA0863	5" Comfort wheel black rubber	<input type="checkbox"/>				<input type="radio"/>
	5" Starec	<input type="checkbox"/>				
DDA0813	techno anthracite	<input type="checkbox"/>	DDA0818	comet red	<input type="checkbox"/>	73,93
DDA0832	artic blue	<input type="checkbox"/>	DDA0819	black	<input type="checkbox"/>	73,93
DDA0862	chrome	<input type="checkbox"/>	DDA0833	copper glaze	<input type="checkbox"/>	73,93
DDA0827	Champagne	<input type="checkbox"/>	DDA0828	Olive	<input type="checkbox"/>	73,93
DDA0817	b-blue	<input type="checkbox"/>				73,93

REAR WHEEL FIXATION

DDA0250	Standard	<input type="checkbox"/>				<input checked="" type="radio"/>
DDA0248	Light rear wheel extension (only positions 1,2,3) (nc camber tube carbon)	<input type="checkbox"/>				180,09
DDA0251	Rear wheel extension positions 1 to 5 (nc nc camber tube carbon and active anti tippers)	<input type="checkbox"/>				228,44

CAMBER AXLE TYPE 1st AXLE

DDA0276	Vario-Ax, adjustable camber 1/3/7/10° (active brakes compulsory with 10° / nc performance brake)	<input type="checkbox"/>				350,71
DDA0277	Standard (aluminium new génération -100g)	<input type="checkbox"/>				<input checked="" type="radio"/>
DDA0278	Carbon = (nc active anti tippers)	<input type="checkbox"/>	-80g			127,01

POSITION OF 1st AXLE

DDA0263	1 (compulsory if using a second axle)	<input type="checkbox"/>				<input type="radio"/>
DDA0264	2	<input type="checkbox"/>				<input type="radio"/>
DDA0265	3	<input type="checkbox"/>				<input checked="" type="radio"/>
DDA0266	4	<input type="checkbox"/>				<input type="radio"/>
DDA0267	5	<input type="checkbox"/>				<input type="radio"/>

CAMBER ANGLE 1st AXLE

DDA0270	0°	<input type="checkbox"/>				<input type="radio"/>
DDA0271	3°	<input type="checkbox"/>				<input checked="" type="radio"/>
DDA0272	6°	<input type="checkbox"/>				<input type="radio"/>

2nd CAMBER AXLE TYPE

DDA0252	Standard	<input type="checkbox"/>				180,09
DDA0280	Vario-Ax	<input type="checkbox"/>				350,71



0338

CPS color set =
Hub
 (wheels 22" and 24" Kuschall moyeu Starec hub and 24" high performance) + **Rim**
 (only for wheels 22" and 24")
+ Forks + Pulgs (Forks/front frame/quick release shaft)
+ Balance adjustable angle backrest

Front wheel not included



0800



0803



0817



0824




0801/0863




0810



0248



Vario-Ax 0276




0251


CAMBER ANGLE 2nd AXLE			
<i>Max. 3° diff. camber with 1st axle - If only one of the axles is a Vario-Axle, 2nd axle only in 0° or 3° camber</i>			
DDA0270	0°	<input type="checkbox"/>	DDA0272 6° <input type="checkbox"/>
DDA0271	3°	<input type="checkbox"/>	<input type="checkbox"/>
QUAD RELEASE AXLE			
DDA0909	Tetra Clip	<input type="checkbox"/>	113,74
REAR WHEEL SIZE (see all the possibilities on Rear wheel type)			
DDA0902	24"	<input type="checkbox"/>	◆
DDA0925	25" (only with wheel High Performance and Spinergy)	<input type="checkbox"/>	○
DDA0903	26" (nc active brakes and active anti tippers)	<input type="checkbox"/>	381,04
REAR WHEEL TYPE (the weight of the wheel pair is indicated)			
DDA0915	Küschall Starec hub (hub CPS if selected)	24", 26" <input type="checkbox"/>	◆
DDA0908	High Performance =(black rim, black hub Starec, hub CPS if selected, crossed spokes, medium pressure tyres or profile ou Marathon Evolution and Speedrun, handrim alu light or titane or supergripp)(nc spoke guard)	24" <input type="checkbox"/>	369,67
DDA0918	High Performance =(black rim, black hub Starec, hub CPS if selected, crossed spokes, tyres Marathon Evolution and Ultrero ZX, handrim alu light)(nc spoke guard) -700g	25" <input type="checkbox"/>	464,45
DDA0906	Spinergy LKX = rim, hub and black spokes) (nc: spokes guard) NEW	24" <input type="checkbox"/>	559,24
DDA0927	High Performance wheels with hub brakes (Starec black hub and silver crossed spokes inside and outside radial) NEW	24" <input type="checkbox"/>	573,46
DDA0913	Mountain Bike (grey rim and Starec hub, crossed spokes, VTT tyres, handrim alu or supergripp)	24" <input type="checkbox"/>	416,11
Rear wheel			
DDA1001	24", 26" Light wheel 7 bar	<input type="checkbox"/>	◆
DDA1002	24" Profile tyre 7 bar	<input type="checkbox"/>	○
DDA1004	24" Solid tyre, profil, grey	<input type="checkbox"/>	○
DDA1006	Bangade trapézoïdal (KIK), noir	<input type="checkbox"/>	56,87
DDA1007	24" Solid tyre black (KIK type)	<input type="checkbox"/>	○
DDA1010	24", 25", 26" "Marathon Evolution" tyre from Schwalbe	<input type="checkbox"/>	56,87
DDA1011	24", 25" Schwalbe speedrun (black and grey) -460g	<input type="checkbox"/>	56,87
DDA1013	24" Schwalbe speedrun (black and blue)	<input type="checkbox"/>	56,87
ANTI PLATT			
DDA1020	Anti-flat tyre tube protection (pair) (nc Marathon Evolution)	<input type="checkbox"/>	76,78
HANDRIM (see all possibilities in "wheel type")			
DDA1113	24", 25" Light Alu dist. 30 mm	<input type="checkbox"/>	◆
DDA1114	24", 25" Light Alu dist. 40 mm	<input type="checkbox"/>	○
DDA1102	24", 26" Supergripp, dist. 30 mm	<input type="checkbox"/>	112,80
DDA1104	24", 26" Steel polished, dist. 30 mm	<input type="checkbox"/>	80,57
DDA1106	24" Titanium, dist. 30 mm	<input type="checkbox"/>	218,01
DDA1120	24", 25" Surge LT, dist. 30 mm	<input type="checkbox"/>	222,75
FRAME OPTION			
SPOKE GUARD (Mountain Bike : 22" only Transparent)			
DDA1200	22" 24", 26" transparent	<input type="checkbox"/>	143,13
DDA1202	24" küschall	<input type="checkbox"/>	143,13



SAFETY			
BRAKES			
DDA1300	Standard brake, lever bent	<input type="checkbox"/>	◆
DDA1302	Performance brake (nc camber 9°)	<input type="checkbox"/>	55,92
DDA1303	Active brake	<input type="checkbox"/>	55,92
DDA1399	Without brakes (only for hub brakes option)	<input type="checkbox"/>	○
SAFETY OPTION			
BRAKE EXTENSION LEVER			
DDA1310	Extension lever for standard brake (Only with levers bent)	<input type="checkbox"/>	55,92
ANTI TIPPERS			
DDA1702	Antitipper tube right	<input type="checkbox"/>	93,84
DDA1703	Antitipper tube left	<input type="checkbox"/>	93,84
DDA1704	Antitipper tube both sides	<input type="checkbox"/>	187,68
DDA1705	Active Antitipper (ground clearance under anti-tip reduced to 10 mm)	<input type="checkbox"/>	218,01
PASSIVE LIGHTS			
DDA1710	Passive lights assembled	<input type="checkbox"/>	16,11
DDA1711	Passive lights enclosed	<input type="checkbox"/>	○
ACCESSORIES			
Backpack accessories			
DDA1744	Package holder	<input type="checkbox"/>	85,31
DDA1745	Bottle holder	<input type="checkbox"/>	36,97
DDA1746	Mobile phone and I phone holder	<input type="checkbox"/>	23,70
DDA1728	Wheel bag 24, 25 et 26"	<input type="checkbox"/>	46,45
TRANSIT WHEELS			
DDA1700	Transit wheels 3" (Only with active antitipper or antitipper tube right)	<input type="checkbox"/>	236,97
Küschall WATCH			
DDA1712	küschall watch, right side	<input type="checkbox"/>	109,00
DDA1713	küschall watch, left side	<input type="checkbox"/>	109,00
TOOLS			
DDA1715	Air pump (nc SB 340mm)	<input type="checkbox"/>	18,01
DDA1716	AirMan	<input type="checkbox"/>	118,48
DDA1717	Toolkit	<input type="checkbox"/>	8,53
FRAME PROTECTORS			
DDA1739	Front frame protectors "küschall" 2pc	<input type="checkbox"/>	33,18
WHEEL BAG including an additionnal wheel pair (without axles)			
DDA1734	24" Mountain bike wheels with bag (with axles)	<input type="checkbox"/>	463,51
SPECIAL BUILT			
DDA0000	Full description and drawing have to be enclosed	<input type="checkbox"/>	



1702 / 1703



1705



1745

1744

1746



1712 / 1713

1728



The Public Prices of options are valid only if they are ordered with the wheelchair

Regarding the options (accessories) which are ordered separately, please consult your Dealer to get the price list

COLLECTION küschall (livrés séparément - remise accessoires)			
1550289	Küshall Sweatshirt black size S	<input type="checkbox"/>	72,50
1550290	Küshall Sweatshirt black size M	<input type="checkbox"/>	72,50
1550291	Küshall Sweatshirt black size L	<input type="checkbox"/>	72,50
1550292	Küshall Sweatshirt black size XL	<input type="checkbox"/>	72,50
1518301	Küshall Sweatshirt grey and black size XXL	<input type="checkbox"/>	72,50
1513357	Küshall Jacket S	<input type="checkbox"/>	145,83
1513358	Küshall Jacket M	<input type="checkbox"/>	145,83
1513359	Küshall Jacket L	<input type="checkbox"/>	145,83
1513360	Küshall Jacket XL	<input type="checkbox"/>	145,83
1513361	Küshall Jacket XXL	<input type="checkbox"/>	145,83
1486614	Küshall umbrella	<input type="checkbox"/>	26,67
1531187	Black cap and silver grey küschall	<input type="checkbox"/>	9,17
1531382	T-Shirt küschall black and silver grey M man	<input type="checkbox"/>	23,33
1531383	T-Shirt küschall black and silver grey L man	<input type="checkbox"/>	23,33
1531463	T-Shirt küschall black and silver grey XL man	<input type="checkbox"/>	23,33
1531464	T-Shirt küschall black and silver grey XXL man	<input type="checkbox"/>	23,33
1531189	T-Shirt küschall black and silver grey S woman	<input type="checkbox"/>	23,33
1531379	T-Shirt küschall black and silver grey M woman	<input type="checkbox"/>	23,33
1531380	T-Shirt küschall black and silver grey L woman	<input type="checkbox"/>	23,33
1531381	T-Shirt küschall black and silver grey XL woman	<input type="checkbox"/>	23,33




Sweat-shirt



Jacket



Pen set



Umbrella



Cap



T-shirt

CROSS COMPATIBILITY

Frame	24" rear wheels			25" rear wheels			26" rear wheels			
	SHv	UL (footrest)		SHv	UL (footrest)		SHv	UL (footrest)		
		Std.	flip-up		Std.	flip-up		Std.	flip-up	
Aluminium	75° Short	45	33-42	29-42	46	33-43	29-43	47	33-44	29-44
		46	34-43	30-43	47	34-44	30-44	48	34-45	30-45
		47	35-44	31-44	48	35-45	31-45	49	35-46	31-46
	75° Standard	48	37-45	33-45	49	37-46	33-46	50	37-47	33-47
		49	38-46	34-46	50	38-47	34-47	51	38-48	34-48
		50	39-47	35-47	51	39-48	35-48	52	39-49	35-49
	75° +5cm	50	39-47	35-47	51	39-48	35-48	52	39-49	35-49
		51	40-48	36-48	52	40-49	36-49	53	40-50	36-50
		52	41-49	37-49	53	41-50	37-50	54	41-51	37-51
	90° Short	45	31-41	27-41	46	31-42	27-42	47	31-43	27-43
		46	32-42	28-42	47	32-43	28-43	48	32-44	28-44
		47	33-43	29-43	48	33-44	29-44	49	33-45	29-45
	90° Standard	48	35-44	31-44	49	35-45	31-45	50	35-46	31-46
		49	36-45	32-45	50	36-46	32-46	51	36-47	32-47
		50	37-46	33-46	51	37-47	33-47	52	37-48	33-48
	90° +5cm	50	37-46	33-46	51	37-47	33-47	52	37-48	33-48
		51	38-47	34-47	52	38-48	34-48	53	38-49	34-49
		52	39-48	35-48	53	39-49	35-49	54	39-50	35-50
	75° Adduction	48	37-45	33-45	49	37-46	33-46	50	37-47	33-47
		49	38-46	34-46	50	38-47	34-47	51	38-48	34-48
		50	39-47	35-47	51	39-48	35-48	52	39-49	35-49
	90° Adduction	48	35-44	31-44	49	35-45	31-45	50	35-46	31-46
		49	36-45	32-45	50	36-46	32-46	51	36-47	32-47
		50	37-46	33-46	51	37-47	33-47	52	37-48	33-48
Titanium	75° Titanium	48	37-45	33-45	49	37-46	33-46	50	37-47	33-47
		49	38-46	34-46	50	38-47	34-47	51	38-48	34-48
		50	39-47	35-47	51	39-48	35-48	52	39-49	35-49
90° Titanium	48	35-44	31-44	49	35-45	31-45	50	35-46	31-46	
	49	36-45	32-45	50	36-46	32-46	51	36-47	32-47	
	50	37-46	33-46	51	37-47	33-47	52	37-48	33-48	
Carbon	Carbon	49	38-46	34-46	50	38-47	34-47	51	38-48	34-48
		50	39-47	35-47	51	39-48	35-48	52	39-49	35-49