

*küschall*<sup>®</sup>  
designed for life



*Go anywhere.  
Do anything.*

*küschall*<sup>®</sup> Ultra-Light



## Technical Data



Seat width:  
28 – 50 cm (in increments of 2 cm)



Backrest angle:  
82°/86°/90°/94°/98°/102°



Seat depth:  
32 – 50 cm (in increments of 2 cm)



Total width  
(with 3° camber):  
SW + 17 cm  
total width folded  
approx. 29 cm



Seat height:  
front: 43 – 53 cm  
rear: 39 – 50 cm  
(in increments of 1 cm)



Total length:  
70°: 88 cm  
80°: 83 cm



Total weight  
unloaded:  
Approx. 9,9 kg



Backrest height:  
30 – 51 cm (in increments of 1.5 cm)



Max. user weight:  
130 kg



Knee-to-heel length:  
20 – 50 cm (in increments of 1 cm)



Transport weight  
w/out rear wheels,  
Approx. 7,5 kg

All measurements and weights indicated are based on seat width 40 cm and seat depth 40 cm in the lightest and most compact configuration possible. Weight and measurements can change with the different configurations.

## Options



Foldable backrest

Front wheel attachment

Frame protection

Mini pushhandles

Adduction frame


Carbon Mudguard

These are only a few of the many options available for this model. Contact your nearest küschall specialist for more information.

## Colours



## CPS (Colour Power Set)



CPS set includes: High Performance wheel hubs with quick release protection, forks, castor housing caps and seat frame caps.



*Good is not good enough.  
It needs to be outstanding.*

**Vanessa Elster,**

*I lead an active lifestyle. Always on the go. My küschall Ultra-Light is practical, light and with a few quick moves it's in the car and I'm on my way.*

*Fold, Lift and Go!*

Küschall AG | Benkenstrasse 260 | CH-4108 Witterswil | Switzerland | Phone +41 (0) 61 487 70 70 | Fax +41 (0) 61 487 70 79 | [www.kuschall.com](http://www.kuschall.com)

**United Kingdom:** Invacare Ltd. | South Road | Bridgend Industrial Estate | Bridgend CF31 3PY | Phone +44 (0)1 656 647327 | Fax +44 (0)1 656 649016 | [www.invacare.co.uk](http://www.invacare.co.uk)

**Ireland:** Invacare Ltd. | Unit 5, Seatown Business Campus | Seatown Road | Swords, Co. Dublin | Phone +353 (0)1 810 7084 | Fax +353 (0)1 810 7085 | [www.invacare.ie](http://www.invacare.ie)

**European Distributor Organisation:** Invacare | Kleiststraße 49 | D-32457 Porta Westfalica | Phone +49 (0)57 31 754 540 | Fax +49 (0)57 31 754 541 | [edo@invacare.com](mailto:edo@invacare.com)

küschall® is a registered trademark. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication, or of parts thereof, as well as any reproduction of images, must have the written consent of Küschall AG.