

## DOLOMITE STEP UP



## DOLOMITE ALPHA



### Models



**Manual**  
Pull out the pin move to required position release the pin.

PRODUCT CODE: 1521667

PRODUCT CODE: 12230-35

### Technical data

Step Up Manual	580/670/760 mm	700/790/880 mm	670 - 1320 mm	340 - 470 mm	620 - 750 mm
Step Up Hydraulic	580/670/760 mm	700/790/880 mm	670 - 1320 mm	340 - 470 mm	620 - 750 mm
Step Up Electric	580/670/760 mm	700/790/880 mm	670 - 1320 mm	340 - 470 mm	620 - 750 mm

  

Step Up Manual	610 mm	16 kg	150 kg	125 mm	215 mm
Step Up Hydraulic	610 mm	18 kg	150 kg	125 mm	215 mm
Step Up Electric	610 mm	23 kg	150 kg	125 mm	215 mm

### Technical data

Basic	650 mm	700 mm	510 mm	950 - 1240 mm	620 mm
Advanced	650 mm	700 mm	510 mm	950 - 1240 mm	620 mm
Basic Platform	650 mm	700 mm	510 mm	1030 - 1240 mm	620 mm
Advanced Platform	650 mm	700 mm	510 mm	1030 - 1240 mm	620 mm

  

Basic	430 mm	200 mm	10.9 kg
Advanced	430 mm	200 mm	11.4 kg
Basic Platform	90 - 270 mm	200 mm	12.4 kg
Advanced Platform	100 - 300 mm	200 mm	12.9 kg

## MANUAL HYDRAULIC PATIENT HOIST

## INVACARE 9805P HYDRAULIC LIFT

## STAND ASSIST POWER HOIST



\* Sling and hanger bar for illustration purposes only



NAME	MANUAL HYDRAULIC PATIENT HOIST	INVACARE 9805P HYDRAULIC LIFT	ELECTRIC – STAND ASSIST HOIST
PRODUCT CODE	APC10100		SMA10601E
SLING	UNIVERSAL	STANDARD SLING	BARIATRIC
FRAME	POWDER COATED MILD STEEL	POWDER COATED FINISH	POWDER COATED MILD STEEL
BASE	63CM – 100CM U BASE IN 5 ADJUSTMENTS	ADJUSTABLE BASE	61CM X 100CM U BASE – ADJUSTABLE
CASTOR	FOUR 4" SOLID CASTORS – REAR WITH BRAKES	5" CASTORS	FOUR 4" SOLID CASTORS – REAR WITH BRAKES
CONTROL	HYDRAULIC HAND PUMP ACTION	HYDRAULIC PUMP	LINAK ACTUATOR OVERRIDE FOR MANUAL RELEASE IN POWER OUT
CARRY CAPACITY	WEIGHT – 150KG GROSS WEIGHT – 47KG	450 LB. GROSS WEIGHT - 63 LB.	LOAD WEIGHT – 175KG NETT WEIGHT – 52KG; GROSS WEIGHT – 60KG
APPLICATION	PREVENT VOCATIONAL INJURY TO NURSE & ENSURE SAFETY OF PATIENT DURING TRANSFER		
BOOM HEIGHT ADJUSTMENT			650 – 2020 MM
CHARGE TIME			4 TO 12 HOURS

**INVACARE  
BIRDIE® FAMILY**



**Birdie**

**Birdie Compact**

Transportability and patient comfort go hand-in-hand. With maximum space for the user, the Invacare Birdie, Birdie Compact and Birdie Plus offer a truly comfortable lift and transfer to or from beds, chairs or even the floor. The lifters are designed to ensure that folding and unfolding can be carried out easily and without the need for tools. In addition, the lifter can be dismantled into two parts if needed – also without tools.



**Manual leg spread**  
Simply step on the pedal to adjust the leg spread.



**Handle for leg spread**  
Handle can easily be attached to operate the manual leg spread.



**Electrical leg spread**  
Available for Invacare Birdie Standard version.



**2 point quick release spreader bars**  
Easy to attach and detach. Available widths: 35 cm, 45 cm, 55 cm.

**Technical data**

<b>Birdie</b>	115 cm	65,5 cm	42 kg	46,5 - 168,5 cm	10 cm
<b>Birdie Compact</b>	100 cm	53 cm	36 kg	51,5 - 158,5 cm	10 cm
<b>Birdie</b>	104 cm	170 kg	40 cycles	2,9 Ah	140 cm
<b>Birdie Compact</b>	89 cm	150 kg	40 cycles	2,9 Ah	107 cm

All measures are taken with 75 mm castors and 2 point spreaderbar. For 100 mm castors add 1,5 cm for height measures and 2 cm for width.

OUTDOOR / POOL HOIST



INVACARE KIRILOS POWER HOIST



Compact and easy to maneuver - the best choice for homecare

Compact in size and offering maximum maneuverability, these mobile lifters provide the ideal solution for the homecare environment. Carers will appreciate the easy handling and smooth operating functions of the Atlante and Kirilos models. The Atlante lift comes in both an electric and hydraulic version; Kirilos is purely electric. Both models come equipped with a handle for easy leg opening.

NAME	OUTDOOR POOL HOIST	INVACARE KIRILOS POWER HOIST
PRODUCT CODE	IGAT180	KIR001
MODEL	IGAT-180	
FRAME	STAINLESS STEEL	MILD STEEL
POWER	WATER PRESSURE	BATTERY 4.5ah
ROTATION	180°	
FOOTREST	FLIP UP FOOTREST INCLUDED	
OPERATION	ALLOWS FOR INDEPENDENT OPERATION	HAND REMOTE
SAFETY	SEAT BELT INCLUDED	
LOAD CAPACITY	160KG WEIGHT – 60KG IN 2 BOXES	140KG