

TOP END OSR RACER & HANDCYCLES



TOP END BASKETBALL & TENNIS



<p>Invacare® Top End® Pro™ Tennis Wheelchair</p>   <p><i>Sleek curved rear bracket offers 6" of quick release rear seat height adjustment.</i></p>	<p>Invacare® Top End® Pro™ Basketball Wheelchair</p>   <p><i>Center-of-gravity and front seat height adjustments are quick and easy.</i></p>	 <p><i>Recessed footrest features toe guard, heel straps, foot strap or optional click strap (shown).</i></p>
--	--	--

TOP END CROSSFIRE ALL TERRAIN



KUCHALL K-SERIES



PRODUCT CODE: CRFAT

USE:	<input checked="" type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input checked="" type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input type="checkbox"/> Folding	<input checked="" type="checkbox"/> Back-height	<input checked="" type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input checked="" type="checkbox"/> COG	<input checked="" type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input checked="" type="checkbox"/> 8 x 2"	<input checked="" type="checkbox"/> S/depth

PRODUCT CODE: KUSKA

USE:	<input checked="" type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input checked="" type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input type="checkbox"/> Folding	<input checked="" type="checkbox"/> Back-height	<input checked="" type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input checked="" type="checkbox"/> COG	<input checked="" type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input checked="" type="checkbox"/> S/depth



Features

- 6061 T6 Aluminium frame
- Adjustable Centre-of-gravity
- Camber 0-9 °
- Lock down folding back
- Adjustable back height
- Adjustable tension seat upholstery
- Adjustable tension back upholstery
- 24" or 25" performance Sun® wheels
- Hard anodised handrims
- 24" or 25" x 2.75" knobby tired
- 8" x 2" front casters
- Adjustable height and angle footplate
- Powder coat paint
- Weight capacity 250 lb.



- 24" or 25" x 2.75" knobby wheels and tires
- 8" x 2" front casters
- Engineered for front and side-to-side stability
- Center-of-gravity adjustment
- Anti-tip option
- Optional push handles, wheel locks and anti-tip bars

Technical Data

Seat width: 240 – 480 mm in increments of 20 mm	Seat depth: 275 – 525 mm in increments of 25 mm	Seat height: front: 450 – 520 mm rear: 380 – 490 mm in increments of 10 mm	Backrest height: 270 – 480 mm in increments of 25 mm
Knee-to-foot length: 290 – 480 cm in increments of 10 mm	Centre of gravity: 80 mm/95 mm 110 mm/125 mm 140 mm	Backrest angle: 24°/28°/32°/38°/50°	Total width: 51W + 170 mm
Total length: 75° approx. 820 mm 90° approx. 760 mm	Total weight, unloaded: approx. 7,7 kg	Max. user weight: 120 kg	Transport weight w/out max wheels: approx. 5,4 kg