

INVACARE AZALEA MINOR



Rea Azalea Minor – designed for petite adults and teenagers
 The Rea Azalea Minor is designed for users that require a seat that has a smaller depth and width. Adapted from the Rea Azalea, the Azalea Minor boasts all the advantages of a reliable, tilt-able wheelchair with its special weight-shifting mechanism. This model is ideal for all users that require a solid, stable wheelchair. Now with new frame colours and upholstery. A tension adjustable backrest is available as a new option.

PRODUCT CODE: BNC0014



Tension adjustable backrest
 A new option adjustable in height and as wide as the seat.



Table tray
 It is transparent and the rail offers gripping opportunities for posture solutions.



Shorter Chair
 The chassis is 30 mm shorter, the seat depth is 50 mm shorter and the backrest 100 mm lower than on standard models. Here in new colour, electric green.

USE:	<input type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input checked="" type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input type="checkbox"/> Medium	<input checked="" type="checkbox"/> Low
ADJUSTABILITY:	<input type="checkbox"/> Folding	<input checked="" type="checkbox"/> Back-height	<input checked="" type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input checked="" type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input checked="" type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

Technical data

Rea Azalea Minor	340 - 440 mm	380 - 450 mm	400/450 mm	550 - 650 mm	240 - 360 mm	330 - 500 mm	+1° to 25°
Rea Azalea Minor	0° - 30°	SW + 220 mm (SW + 230 mm)	900 - 1250 mm	900 - 970 mm	32 kg	Max 75 kg	Transport weight 14,5 kg
Rea Azalea Minor	Crash-tested ISO 7176-19						

CPC WITH TILT AND RECLINE

TRANSPORTER BUGGY



Model ALK958LC-36

- 1) Aluminium lightweight frame with comfortable seat, reclining high back, adjustable headrest
- 2) 6" PVC castor 16" wheel
- 3) Adjustable and detachable armrest, elevating footrest, seat and backrest adjustable
- 4) With elasticity belt

Weight – 12kg
 Load Capacity – 50kg
 Sizes – 12", 14", 16"
 Frame – Mild Steel
 Rear Tyre – 11" Solid

USE:	<input type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input checked="" type="checkbox"/> Passive
ENVIRONMENT:	<input checked="" type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input type="checkbox"/> Medium	<input checked="" type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input checked="" type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input checked="" type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

USE:	<input type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input checked="" type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input type="checkbox"/> Medium	<input checked="" type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input checked="" type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

NAME	CPC WITH TILT AND RECLINE	TRANSPORTER BUGGY
PRODUCT CODE	ALK958LC / YK9140	BUG001
SIZE	12", 14", 16", 18"	12", 14", 16"
FRAME	ALUMINIUM LIGHTWEIGHT	MILD STEEL
TYRE	6" PVC CASTOR 16" WHEEL	11" SOLID REAR WHEEL
LOAD CAPACITY	50KG (18" - 60KG)	50KG
WEIGHT	15KG	12KG

PAEDIATRIC REBOTEC WALKING AIDS



PRODUCT CODE: 1866020

FOX

Aluminium rollator for children/youth/woman

Equipment:

- dual action brake
- foldable
- airless tyres
- particular light through powder coated aluminium frame

ref. no. 186...20



PRODUCT CODE: 1870010/20



PLUTO

Children reverse walker

- the PLUTO reverse walker provides independent movement for children.
- it supports the user from behind and promotes upright posture when standing and walking.
- it supports the therapeutic exercises to learn natural movements and increase the balance.
- the braking system in the back wheels prevents unintentional rolling back.



WALK-ON

Aluminium walking frame

Equipment:

- particular light through powder coated aluminium frame
- nonslipping handles for safer use
- easy foldable

ref. no. 174...



PRODUCT CODE: 2872010

PRODUCT CODE: 1740

ref. no. 174.10



WALK-ON

Aluminium walking frame with 2 rolls

Equipment:

- particular light through powder coated aluminium frame
- nonslipping handles for safer use
- easy foldable

ref. no. 175...

ref. no. 175.11

PRODUCT CODE: 17420


SAFE-IN-KIDS

Crutch with nearly closed forearm cuff



PRODUCT CODE: 1229510

PAEDIATRIC ONE TOUCH WALKER



- Anodized Aluminium Frame – Silver
- 1" Diameter Tubing
- Black crossbrace and accessories
- Quantity – 4 per container

DUAL RELEASE PADDLE WALKER



REVERSE WALKER

ALK763L Walker

					N.W	PCS	G.W
57	-	73	47-65	360x25	4	2	10



NAME	ONE TOUCH WALKER	NAME	DUAL RELEASE PADDLE WALKER
PRODUCT CODE	12852S	PRODUCT CODE	6291-JR
ADJUSTABLE HEIGHT	28" TO 32"	WEIGHT	LOAD CAPACITY – 300lb. GROSS WEIGHT – 5lb1oz
FRAME	1" DIAMETER TUBING	PATIENT HEIGHT	4'4" - 5'7"
COLOUR	BLACK CROSSBRACE AND ACCESSORIES	HEIGHT ADJ.	25.375" - 32.375"
PACKING	4 PER CONTAINER	WIDTH INSIDE HANDGRIPS	17"
WIDTH INSIDE HANDGRIPS	17"	LIMITED WARRANTY	5 YEARS

KIDDIES SATURN POWER



KIDDIES OMEGA POWER



USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input checked="" type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input checked="" type="checkbox"/> H/Adj arm	<input checked="" type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input checked="" type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input checked="" type="checkbox"/> H/Adj arm	<input checked="" type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

NAME	SATURN KIDDIES POWER	OMEGA A1 POWER
PRODUCT CODE	SATPOWK	A1POWK
SEAT	12" AND 14"	14"
CASTORS	8X2" SOLID	8X2" SOLID
REAR TYRE	12½ BY 2¼ " SOLID	12½ X 2¼ " FOAM FILLED
CONTROL	DYNAMIC REM D01	DYNAMIC REM D01
BATTERY	12V36A/H	12V36A/H
ARMREST AND FOOT	REMOVABLE & HEIGHT ADJUSTABLE ARM REMOVABLE FOOTREST	REMOVABLE & HEIGHT ADJUSTABLE ARM REMOVABLE FOOTREST

BASIK KIDDIES

SATURN KIDDIES



USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input checked="" type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input checked="" type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

NAME	BASIK KIDDIES	SATURN KIDDIES
PRODUCT CODE	YK9031K	KIDDIES
SIZE	84CM	12"
COLOURS	MULTICOLOURED FRAME	WET BLACK
FRAME	FOLDING MILD STEEL POWDER COATED	FOLDING MILD STEEL POWDER COATED
ARMS	DESK, FLIP BACK	HEIGHT ADJUSTABLE FULL ARM
FOOTREST	REMOVABLE	REMOVABLE 90° FOOTREST
CASTORS	6" CASTOR	8X2" PU CASTOR ON STEEL FORK
RIMS	20" REAR WHEEL	20 X 1.75 STEEL SPOKE RIM WITH 24X1.75 TYRE WITH NO MORE FLAT (NMF), FIXED AXLES
UPHOLSTERY	BLACK NYLON	BLACK REINFORCED NYLON
WEIGHT	LOAD CAPACITY - 60KG GROSS WEIGHT - 18.5KG	LOAD CAPACITY - 50KG GROSS WEIGHT - 16KG
PACKING	80X34X85CM	100X34X85CM

SATURN ALU-LITE KIDDIES



INVACARE ACTION 3 JUNIOR



USE:	<input checked="" type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input checked="" type="checkbox"/> COG	<input checked="" type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input checked="" type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

USE:	<input checked="" type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input checked="" type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input checked="" type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

NAME	SATURN ALU-LIGHT	INVACARE ACTION 3 JUNIOR
PRODUCT CODE	LW25-77K	CNB0030
SIZE	14"	CHILD 8"-1" JUNIOR 12"-15"
COLOURS	METALLIC BLUE, WET BLACK	OPTIONAL
FRAME	FOLDING ALUMINIUM POWDER COATED	FOLDING ALUMINIUM FRAME
ARMS	FIXED HEIGHT FLIP BACK FULL ARM	REMOVABLE & HEIGHT ADJUSTABLE "T" ARM
FOOTREST	REMOVABLE	REMOVABLE
CASTORS	8X1" SOLID CASTOR ON ALUMINIUM FORK	6" SOLID CASTER
RIMS	24X13/8 MAG RIM WITH PU TYRE 12MM QUICK RELEASE AXLES	20" ALUMINIUM SPOKE RIM WITH PU TYRE QUICK RELEASE AXLES
UPHOLSTERY	BLACK REINFORCED NYLON	BLACK REINFORCED NYLON
ADJUSTMENT	CENTRE OF GRAVITY, LOW SLUNG HEIGHT ADJUSTABLE BACK REST WITH PUSH HANDLES	CENTRE OF GRAVITY LOW SLUNG
WEIGHT	LOAD CAPACITY – 115KG GROSS WEIGHT – 17KG (100X34X85CM)	LOAD CAPACITY – 80KG GROSS WEIGHT – 15KG

KIDDIES GROWING WHEELCHAIR



Seat Additional Strap



Rear Axle Close Up

USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input checked="" type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input checked="" type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input checked="" type="checkbox"/> Seat angle	<input checked="" type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input checked="" type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input checked="" type="checkbox"/> E/leg rest
	<input checked="" type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

NAME	KIDDIES GROWING WHEELCHAIR
PRODUCT CODE	KIDDIESGR
SEAT WIDTH	10", 12", 14"
SEAT DEPTH	10", 12", 14"
COLOURS	BLACK
FRAME	POWDER COATED ALUMINIUM
CASTOR	6" SOLID
RIMS	20" ALUMINIUM SPOKE RIM
ADJUSTMENT	CENTRE OF GRAVITY, SEAT HEIGHT
UPHOLSTERY	BLACK NYLON
WEIGHT	12KG

CTM MICRO SCOOTER



USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input checked="" type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input checked="" type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

CTM MICRO SCOOTER – HS118

Type	Mobility Scooters	Curb Climbing	35mm/1.4"
Model Number	HS-118	Turning Radius	1270mm/50"
Overall Length	1000mm/39.4"	Suspension	N/A
Overall Width	510mm/20"	Brake	Electro-Mechanical
Overall Height	890mm/35"	Seat Type	Padded Foldable
Wheels: Front	195mm/8"	Seat Width	425mm/17"
Wheels: Rear	195mm/8"	Motor Size	250W, 4700 rpm
Weight w/ Batteries	40kg/88lbs	Battery Size	(2) 12V • 12Ah
Max. Speed	6.4kmph/4mph	Battery Weight	9.2kg/20lbs
Weight Capacity	115kg/253lbs	Travel Range	10km/6.2miles
Ground Clearance	30mm/1.2"	Battery Charger	1.8A Off Board
Grade Climbable	8 degrees	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Blue		

CTM MICRO SCOOTER – HS125

Type	Mobility Scooters	Curb Climbing	40mm/1.6"
Model Number	HS-125	Turning Radius	10055mm/9.6"
Overall Length	1000mm/39.4"	Suspension	N/A
Overall Width	470mm/18.5"	Brake	Electro-Mechanical
Overall Height	880mm/34.6"	Seat Type	Padded Foldable
Wheels: Front	190mm/8"	Seat Width	425mm/16.7"
Wheels: Rear	190mm/8"	Motor Size	250W 4600 rpm
Weight w/ Batteries	37kg/81.6lbs	Battery Size	(2) 12V • 12Ah
Max. Speed	8kmph/5mph	Battery Weight	9.3kg/20.35lbs
Weight Capacity	115kg/253lbs	Travel Range	9km/5.6 miles
Ground Clearance	35mm/1.4"	Battery Charger	2A On Board
Grade Climbable	8 degrees	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Blue		

CTM MINI SCOOTER



USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input checked="" type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

USE:	<input type="checkbox"/> Active	<input checked="" type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input checked="" type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input checked="" type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

CTM MINI SCOOTER – HS295

Type	Mobility Scooters	Curb Climbing	30mm/1.2"
Model Number	HS-295	Turning Radius	1270mm/50"
Overall Length	990mm/39"	Suspension	N/A
Overall Width	505mm/19.9"	Brake	Electro-Mechanical
Overall Height	870mm/34.3"	Seat Type	Padded Foldable Swivel
Wheels: Front	190mm/8"	Seat Width	400mm/15.8"
Wheels: Rear	190mm/8"	Motor Size	250W 4600 rpm
Weight w/ Batteries	46.76kg/103lbs	Battery Size	(2) 12V • 12 Ah
Max. Speed	8kmph/5mph	Battery Weight	10.4kg/23lbs
Weight Capacity	115kg/250lbs	Travel Range	10km/6.2 miles
Ground Clearance	40mm/1.6"	Battery Charger	2A On Board
Grade Climbable	8 degrees	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Silver, Blue		

CTM MINI SCOOTER – HS265

Type	Mobility Scooters	Curb Climbing	50mm/2"
Model Number	HS-265	Turning Radius	1060mm/41"
Overall Length	1010mm/39.8"	Suspension	N/A
Overall Width	505mm/19.9"	Brake	Electro-Mechanical
Overall Height	900mm/35.4"	Seat Type	Padded Foldable Swivel
Wheels: Front	190mm/8"	Seat Width	430mm/17"
Wheels: Rear	190mm/8"	Motor Size	250W 4600 rpm
Weight w/ Batteries	46kg/101.4lbs	Battery Size	(2) 12V • 12 Ah
Max. Speed	8kmph/5mph	Battery Weight	10.4kg/23lbs
Weight Capacity	115kg/250lbs	Travel Range	10km/6.2 miles
Ground Clearance	40mm/1.6"	Battery Charger	2A On Board
Grade Climbable	8 degrees	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Silver, Blue, Turquoise		

CTM MID-LINE SCOOTER



USE:	<input checked="" type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input checked="" type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

USE:	<input checked="" type="checkbox"/> Active	<input type="checkbox"/> Med-Active	<input type="checkbox"/> Passive
ENVIRONMENT:	<input type="checkbox"/> Indoor	<input type="checkbox"/> Indoor/Outdoor	<input checked="" type="checkbox"/> Outdoor/Indoor
ACTIVITY LEVEL:	<input type="checkbox"/> High	<input checked="" type="checkbox"/> Medium	<input type="checkbox"/> Low
ADJUSTABILITY:	<input type="checkbox"/> Folding	<input type="checkbox"/> Back-height	<input type="checkbox"/> Back angle
	<input type="checkbox"/> Seat angle	<input type="checkbox"/> COG	<input type="checkbox"/> Q/R wheel
CONFIGURABILITY:	<input type="checkbox"/> TAB	<input type="checkbox"/> Full arm	<input type="checkbox"/> E/leg rest
	<input type="checkbox"/> H/Adj arm	<input type="checkbox"/> 8 x 2"	<input type="checkbox"/> S/depth

CTM MID-LINE SCOOTER – HS320

Type	Mobility Scooters	Curb Climbing	50mm/2"
Model Number	HS-320	Turning Radius	1020mm/40.2"
Overall Length	1090mm/42.9"	Suspension	N/A
Overall Width	550mm/21.7"	Brake	Electro-Mechanical
Overall Height	870mm/34.3"	Seat Type	Light Weight Foldable Swivel
Wheels: Front	230mm/9"	Seat Width	460mm/18"
Wheels: Rear	230mm/9"	Motor Size	450W 3700 rpm
Weight w/ Batteries	53.8kg/118.7lbs	Battery Size	(2) 12V • 18 Ah
Max. Speed	8kmph/5mph	Battery Weight	12.7kg/28lbs
Weight Capacity	135kg/300lbs	Travel Range	135km/8.4 miles
Ground Clearance	40mm/1.6"	Battery Charger	2A Off Board
Grade Climbable	12 degrees	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Silver, Blue, Turquoise		

CTM MID-LINE SCOOTER – HS360

Type	Mobility Scooters	Curb Climbing	50mm/2"
Model Number	HS-360	Turning Radius	1415mm/55.7"
Overall Length	1160mm/45.7"	Suspension	N/A
Overall Width	550mm/21.7"	Brake	Electro-Mechanical
Overall Height	870mm/34.3"	Seat Type	Light Weight Foldable Swivel
Wheels: Front	230mm/9"	Seat Width	460mm/18"
Wheels: Rear	230mm/9"	Motor Size	450W 3700 rpm
Weight w/ Batteries	67kg/148.9lbs	Battery Size	(2) 12V • 18Ah
Max. Speed	8kmph/5mph	Battery Weight	12.7kg/28lbs
Weight Capacity	135kg/300lbs	Travel Range	13km/8.1 miles
Ground Clearance	40mm/1.6"	Battery Charger	2A Off Board
Grade Climbable	12 degrees	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Silver, Blue, Turquoise		