

**SINGLE FOLD ALUMINIUM STRETCHER**



**DOUBLE FOLD STRETCHER**



**SCOOP STRETCHER**



**CHAIR /EVACUATION STRETCHER**



| NAME          | SINGLE FOLD STRETCHER            | DOUBLE FOLD STRETCHER            | SCOOP STRETCHER   | CHAIR STRETCHER     |
|---------------|----------------------------------|----------------------------------|-------------------|---------------------|
| PRODUCT CODE  | JL-1                             | JL-2                             | JL-4              | JL-5                |
| FRAME         | ALUMINIUM ALLOY                  | ALUMINIUM ALLOY                  | ALUMINIUM ALLOY   | ALUMINIUM ALLOY     |
| UNFOLDED      | 185LX50WX14H (CM)                | 208LX54.5WX13H (CM)              |                   | 90LX50WX91H (CM)    |
| FOLDED        | 92.5LX50WX9H (CM)                | 104LX9W17H (CM)                  |                   | 90LX59WX17H (CM)    |
| UPHOLSTERY    | BLUE CANVAS                      | BLUE CANVAS                      |                   | ORANGE PVC          |
| GROSS WEIGHT  | 14KG FOR 2 UNITS                 | 14KG FOR 2 UNITS                 | 11KG              | 8KG                 |
| LOAD CAPACITY | 125KG                            | 125KG                            | 140KG             | 140KG               |
| PACKING       | 98.5LX53WX19H (CM)<br>2 PER PACK | 9108LX22WX19H (CM)<br>2 PER PACK | 122Lx45Wx11H (CM) | 93LX54WX17H (CM)    |
| MAX SIZE      |                                  |                                  | 202LX42WX6H (CM)  |                     |
| MIN SIZE      |                                  |                                  | 168LX42WX6H (CM)  |                     |
| ACCESSORIES   |                                  |                                  |                   | CARRY / HANGING BAG |

**JL-9 SELF LOADING  
AMBULANCE STRETCHER**




|                      |                                    |
|----------------------|------------------------------------|
| <b>NAME</b>          | AMBULANCE STRETCHER – SELF LOADING |
| <b>PRODUCT CODE</b>  | JL-9                               |
| <b>FRAME</b>         | ALUMINIUM                          |
| <b>GROSS WEIGHT</b>  | GROSS WEIGHT – 45KG                |
| <b>LOAD CAPACITY</b> | 150KG                              |
| <b>ACCESSORIES</b>   | PVC COVERED FOAM MATTRESS          |
| <b>HIGH POSITION</b> | 1930LX550WX850H (MM)               |
| <b>LOW POSITION</b>  | 1900LX600WX240H (MM)               |
| <b>BACK INCLINE</b>  | UP TO 60°                          |
| <b>CASTORS</b>       | 4XSWIVEL CASTOR                    |

**3 FUNCTION ELECTRICAL HI-LO HOSPITAL BED**

**Features:**

- 1) Cold steel plate wholly moulded surface
- 2) Integrated stamping frame and treated with electrostatic spraying, anti-aging and anti-rust
- 3) ABS bed head, 6 aluminium bars fence can go up and down
- 4) 3pcs electro-motor, 1pc control box
- 5) One set of handset
- 6) 5" universal castors (2pcs with brake, 2pcs without brake)
- 7) Max Loading: 280kg
- 8) Battery backup
- 9) Has been assessed and conforms to the standards of ISO13485:2003



**Dimensions:**

- 1) Size: L2130 x W980 x H450 – 700mm
- 2) Angle of back section: 0 ~ 80° (± 5°)
- 3) Angle of leg section: 0 ~ 40° (± 5°)
- 4) Adjustable height: 450 – 700mm

**Accessories:**

- 1) 1 IV-Pole
- 2) Wheels
- 3) Aluminium alloy protective railing



PRODUCT CODE: ALK06-B02P



3pcs electro motor, 1pc control box one set of handset, battery backup

Can change the plug and voltage of motor as required



Two Sections IV Pole (Four Hooks) Frame: Stainless Steel

5" universal castors (2pcs with brake, 2pcs without brake)

INVACARE MEDLEY ERGO



**Shipping and Storage**  
Medley Ergo beds are very easy to store and transport when placed on the dedicated transporter.

PRODUCT CODE: ERGO


Technical data

|                        |  |   |  |   |                                     |                        |
|------------------------|--|---|--|---|-------------------------------------|------------------------|
|                        |  |   |  |   |                                     |                        |
| <b>Medley Ergo</b>     | Width outside: 1020 mm<br>Width inside: 900 mm | Length outside: 2240 mm<br>Length inside: 2060 mm | Height adjustment †: 400 - 600 mm/<br>330 - 730 mm | Mattress support dimensions: 800-230-300-670 mm | 0 - 70°<br>Backrest angle           | 0 - 24°<br>Thigh angle |
| <b>Medley Ergo Low</b> | Width outside: 1020 mm<br>Width inside: 900 mm | Length outside: 2240 mm<br>Length inside: 1060 mm | Height adjustment †: 290 - 690 mm/<br>210 - 610 mm | Mattress support dimensions: 800-230-300-670 mm | 0 - 70°<br>Backrest angle           | 0 - 24°<br>Thigh angle |
|                        |  |   |  |   |                                     |                        |
| <b>Medley Ergo</b>     | 0 - 11°<br>Legrest lift                        | 11°<br>Tilt option ‡                              | 64 kg  | SWL 180 kg<br>Max. user weight 145 kg           | Heaviest part of the product: 17 kg |                        |
| <b>Medley Ergo Low</b> | 0 - 11°<br>Legrest lift                        | 11°<br>Tilt option ‡                              | 64 kg  | SWL 180 kg<br>Max. user weight 145 kg           | Heaviest part of the product: 17 kg |                        |

**3 FUNCTION MANUAL HI-LO HOSPITAL BED**

**Features:**

- 1) 3sets manual crank system
- 2) ABS bed head, 6 aluminium bars fence can go up and down
- 3) 4pcs soft crashworthy bumpers
- 4) Accessories: 5" universal castors (2pcs with brake, 2pcs without brake), 2pcs urine hooks, 6pcs IV drip holes
- 5) Integrated stamping frame and treated with electrostatic spraying, anti-aging and anti-rust
- 6) Max Loading: 280kg
- 7) Has been assessed and conforms to the standards of ISO13485:2003



**Dimensions:**

- 1) Size: L2130 x W980 x H450 – 700mm
- 2) Angle of back section: 0 - 80° (± 5°)
- 3) Angle of leg section: 0 - 40° (± 5°)
- 4) Adjustable height: 450 – 700mm

**Accessories:**

- 1) 1 IV-Pole
- 2) Wheels
- 3) Aluminium alloy protective railing



PRODUCT CODE: ALK06-A328P



Triple manual crank system  
The crank adapt Hidden design




Two Sections IV Pole (Four Hooks)  
Frame: Stainless Steel  
5" universal castors (2pcs with brake, 2pcs without brake)

**5 FUNCTION ELECTRIC  
TRENDELENBERG HOSPITAL BED**



**Features:**

- 1) Cold steel plate wholly moulded surface, connected with arc soft joint
- 2) Streamline bed headboard and banister made of PE
- 3) 4pcs electro-motor, 1pc control box, 1pc handset
- 4) Automatic lifting railing, safe and reliable, and can be fixed upward and downward
- 5) One set of handset
- 6) Four wheels, two with brakes, two without brakes
- 7) Battery backup
- 7) Has been assessed and conforms to the standards of ISO13485:2003, CE



**Dimensions:**

- 1) Size: L2240 x W1100 x H450-720mm
- 2) Angle of back section: 0 - 80° (± 5°)
- 3) Angle of leg section: 0 ~ 40° (± 5°)
- 4) Angle of Trendelenburg and reserve Trendelenburg: 0 ~ 13° (± 2°)
- 5) Adjustable height: 450 – 720mm

**Accessories:**

- 1) 1 IV-Pole
- 2) Wheels
- 3) Automatic lifting railing



PRODUCT CODE: ALK06-B01A



4pcs electro-motor, 1pc control box, one set of handset, Battery backup

Can change the plug and voltage of the motor as inquiry



Two Sections IV Pole (Four Hooks)  
Frame: Stainless Steel

5" universal castors (2pcs with brake, 2pcs without brake)

**1 FUNCTION MANUAL HOSPITAL BED**


**Features:**

- 1) Single manual crank system
- 2) ABS bed head, 6 aluminium bars fence can go up and down
- 3) 4pcs soft crashworthy bumpers
- 4) Accessories: 5" universal castors (2pcs with brake, 2pcs without brake), 2pcs urine hooks, 6pcs IV drip holes
- 5) Integrated stamping frame and treated with electrostatic spraying, anti-aging and anti-rust
- 6) Max Loading: 280kg
- 7) Has been assessed and conforms to the standards of ISO13485:2003



**Dimensions:**

- 1) Size: L2130 x W980 x H500mm
- 2) Angle of back section: 0 - 80° (± 5°)



**Accessories:**

- 1) 1 IV-Pole
- 2) Wheels
- 3) Aluminium alloy protective railing
- 4) Oddment hook



PRODUCT CODE: ALK06-A132P



Single manual crank system  
The crank adapt Hidden design




Two Sections IV Pole (Four Hooks)  
Frame: Stainless Steel  
5" universal castors (2pcs with brake, 2pcs without brake)

**3 FUNCTION SEMI-ELECTRIC HOSPITAL BED**

**Features:**

- 1) Three functions of the Semi-electric hospital bed, including head and foot positioning and bed height
- 2) This semi-electric hospital bed combines electric motor and the manual. Pendant control and Manual crank for this bed provides raises and lowers the upper body, knees, bed-frame height separately
- 3) Cold steel plate wholly moulded surface. Integrated stamping frame and treated with electrostatic spraying, anti-aging and anti-rust.
- 4) ABS bed head, 6 aluminium bars fence can go up and down
- 5) 3pcs electro-motor, 1pcs control box, one set of handset and crank set
- 6) 5" universal castors (2pcs with brake, 2pcs without brake)
- 7) Max loading 280kg
- 8) Battery backup (option)
- 9) Has been assessed and conforms to the standards of ISO13485-2003, CE



**Dimensions:**

- 1) Size: L2130 x W980 x H450-680mm
- 2) Angle of back section: 0 - 80° (± 5°)
- 3) Angle of leg section: 0 - 40° (± 5°)
- 4) Adjustable height: 450 – 680mm

**Accessories:**

- 1) 1 IV-Pole
- 2) Wheels
- 3) Aluminium alloy protective railing

PRODUCT CODE: ALK06-B07P



This semi-electric hospital bed combines electric motor and the manual. Pendant control and Manual crank for this bed.



Two Sections I.V. Pole (Four Hooks)  
Frame: Stainless steel  
5" universal castors (2pcs with brake, 2pcs without brake)



ALUMINIUM COTSIDE



PRODUCT CODE: ALK06-P302P

TRAUMA TROLLEY



PRODUCT CODE: TRMTRL

EXAMINATION COUCH



PRODUCT CODE: ALK08-MB01  
Aluminium banister

Specifications:  
Overall Dimensions: L1800mm x W600mm x H650mm  
Max. Weight Capacity: 170kgs  
Mattress: High Class Sponge Artificial Leather  
Frame: Stainless Steel

STAINLESS STEEL STRETCHER



PRODUCT CODE: ALK06-A101B  
Stainless Steel Stretcher

Material: stainless steel  
Stainless steel rail can up and down  
Bed frame can be separated  
Can be made as assembling  
Size: L1935 x W570 x H650mm

**HYDRAULIC STRETCHER**

PRODUCT CODE: ALK06-H802



**Features:**

- 1) Size: 2100 x 690 x 460/800mm
- 2) The whole stretcher is epoxy painted with wire surface
- 3) Backrest adjusted by manual crank system, imported hydraulic pump adjust height between 460-850mm
- 4) Strong and useful side rail – four bumpers
- 5) The surface is made of special material which can be used for X-ray
- 6) Bedrail, total length 1450mm, aluminium alloy hand bar with oxidation treatment, six stainless side bars consisted in the side rail
- 7) One pc imported hydraulic pump

**Functions:**

- 1) Angle of back section: 0 - 80° (± 5° )
- 2) Foot rest Trend: 0 - 11° (± 2° )
- 3) Adjusted height: 460 – 800mm

**Accessories:**

- 1) 1 V Pole
- 2) Wheels
- 3) Automatic aluminium railing
- 4) 1 pc imported hydraulic pump
- 5) And push handle

**PATIENT TROLLEY**

PRODUCT CODE: ALK06-A103B



**Features:**

Size: 1950 x 570 x 750cm

- 1) 2-section: adjustable head section and fixed trunk section with mattress
- 2) Head section is adjusted by manual rack mechanism
- 3) Frame is made of tubular stainless steel
- 4) Fitted with 4 anti-shock resin bumpers, one at each corner
- 5) Water proof mattress cover, high dense sponge inside
- 6) Fitted with 2 foldable stainless steel side rails, lockable in raised position
- 7) With multiple drainage bag hooks at each side
- 8) IV pole with clamp, to be mounted at head end
- 9) Fitted with 1 storage basket and 1 oxygen cylinder carrier
- 10) Mounted on 4 anti-static swivel castors, diameter 12.5cm with brakes

**CHILDREN'S COT BED**



PRODUCT CODE: ALK06-CHIP



**Features:**

- 1) Single manual crank system
- 2) Frame is made of mild steel
- 3) Totally enclosed bedside rail
- 4) Bedside rail is possible to be adjusted to three altitudes

**Dimensions:**

- 1) Size: L500 x 800 x 550mm
- 2) Angle of back section: 0 - 80° (± 5°)

**BABY CABINET**

PRODUCT CODE: ALK06-BABY



Baby cabinet in stainless steel material and bassinette in Acrylic plastic  
 Size: L820 x W500 x H950mm  
 Can be made as CKD

**BEDSIDE LOCKER**

PRODUCT CODE: ALK06-AG01



**Features:**

- 1) 2pcs towel stands
- 2) 4pcs oddment hooks
- 3) 1pc water bottle stand
- 4) Optional choose 4pcs 2 castors
- 5) Size: 480 ##### (cut off)
- 6) All colours ##### (cut off)
- 7) Can be ##### (cut off)
- 8) Material ##### (cut off)

**STAINLESS STEEL NURSING CART**

PRODUCT CODE: ALK07-H04



**Structure and configuration:**

- 1) All stainless steel structure, fold-deck with silent wheels and dirt bags.
- 2) Can be made as CKD

**Size:**

L900 x W500 x H900/970mm

**STAINLESS STEEL KICK BUCKET**

PRODUCT CODE: ALK06-H06



Size: 490 x 490 x 380mm

- 1) Stainless steel bucket, capacity 13 litres
- 2) Mounted on a minimum of 4 swivel castors, diameter 5cms

**BEDSIDE TABLE**

PRODUCT CODE: ALK06-AT3



Movable, can adjust the position low or high  
 Made of ABS and aluminium alloy  
 Specifications: L800 X W400 X H750 – 1050MM

**4 SECTION MEDICAL SCREEN**

PRODUCT CODE: ALK06-501



Material: stainless steel  
 Per section: 500mm 1800mm

**ONE STEP FOOTSTOOL**

PRODUCT CODE: ALK06-D01



**Description:**

- 1) Made of steel with powder coated, nonslip tread and mounted on rubber feet

Size:  
 490 X 290 X 210mm