

WALKING TUTOR - HEAVY DUTY

- New aluminium oval tube frame construction.
- For institutions, rental or domestic use.
- Safe working load of 180kg, Australian standards certified.
- Fitted with new "light touch" hand brakes with a lock brake mechanism.
- The frame is widened between the rear wheels for maximum walking space.
- The Gutter arm walker is foldable for transportation.
- The oval tube provides a directional height adjustment function (Press button adjustment).
- The handles bars are multi adjustable for angle and fore/aft function.

**Registered design
9210A**



W/Lock brake Swivel castors

9210A

9210B

overall height	overall length	overall width	armrest length	folded width	castor size	unit weight	weight capacity
min103cm max127cm	88cm	65cm	min30cm max43cm	28cm	20cm	12.2kg	180kg

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za