

ALPHA A6 POWER WITH TILT



NAME	ALPHA TILT POWER
PRODUCT CODE	A6POW
SIZE	14", 16", 18"
BATTERIES	2 X 36AH/12 VOLT
MOTORS	250 WATT
CASTORS	8 X 2" FOAM FILLED GREY
SPEED	6KM/H
RANGE	20KM (TERRAIN AND USER DEPENDENT)
REAR RIMS & TYRES	12 1/2 X 2 1/4 ALUMINIUM RIM 12 1/2 X 2 1/4 ALUMINIUM RIM
ARMS	HEIGHT ADJUSTABLE T ARM
GRADIENT	12 DEGREES
CONTROLS	DYNAMIC SHARK REMD11 POWER MODULE PMB11
SEAT FUNCTION	
LOAD CAPACITY	LOAD CAPACITY – 120KG WEIGHT – 70KG
PACKING	70X70X100CM

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za