

BASIK HYBRID GRANNY



- Wheelchair/Transport chair two way use
- Aluminium chair frame
- Flip back desk arm
- Detachable mag wheel with PU tyre
- 8" rear wheel & 8" front castor
- Tool-free adjustable swing away footrest
- Aluminium front fork with 2 adjustable positions

NAME	BASIK HYBRID GRANNY
PRODUCT CODE	Yk9061
SIZE	18"
COLOURS	SILVER
FRAME	ALUMINIUM CHAIR FRAME
ARMS	FLIP BACK DESK ARM
FOOTREST	TOOL-FREE ADJUSTABLE SWING AWAY FOOTREST
CASTORS	24" REAR WHEEL & 8" FRONT CASTOR
RIMS	FIXED MAG WHEEL WITH PU TYRE
UPHOLSTERY	BLACK NYLON
WEIGHT	LOAD CAPACITY 90KG & WEIGHT 16KG
PACKING	100 x 34 x 85 (cm)
USE	LIGHT INDOOR / SMOOTH SURFACE USE

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za