

## BASIK KIDDIES



<b>NAME</b>	BASIK KIDDIES
<b>PRODUCT CODE</b>	YK9031K
<b>SIZE</b>	12", 14"
<b>COLOURS</b>	MULTICOLOURED FRAME
<b>FRAME</b>	FOLDING MILD STEEL POWDER COATED
<b>ARMS</b>	DESK, FLIP BACK
<b>FOOTREST</b>	REMOVABLE
<b>CASTORS</b>	6" CASTOR
<b>RIMS</b>	20" REAR WHEEL
<b>UPHOLSTERY</b>	BLACK NYLON
<b>WEIGHT</b>	LOAD CAPACITY - 60KG GROSS WEIGHT - 18.5KG
<b>PACKING</b>	80X34X85CM

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za