

## BASIK TRANSIT



<b>NAME</b>	BASIK TRANSIT WHEELCHAIR
<b>PRODUCT CODE</b>	YK9092/3
<b>SIZE</b>	17" & 19"
<b>COLOURS</b>	BLACK
<b>FRAME</b>	MILD STEEL & ALUMINIUM POWDER COATED
<b>ARMS</b>	FIXED FULL ARMS
<b>FOOTREST</b>	REMOVABLE
<b>CASTORS</b>	FOUR 8X1" SOLID CASTOR
<b>RIMS</b>	
<b>UPHOLSTERY</b>	BLACK NYLON
<b>WEIGHT</b>	LOAD CAPACITY – 90KG GROSS WEIGHT – 11KG
<b>PACKING</b>	76 X 26 X 60 (CM)
<b>USE</b>	LIGHT INDOOR / SMOOTH SURFACE USE

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za