

BASIK WHEELCHAIR



NAME	BASIK WHEELCHAIR
PRODUCT CODE	YK9031
SIZE	16", 18", 20"
COLOURS	BLACK
FRAME	MILD STEEL POWDER COATED
ARMS	FLIP BACK DESK ARMS
FOOTREST	REMOVABLE WITH PLASTIC FOOTPLATE
CASTORS	8X1" SOLID CASTOR
RIMS	24X1" SPOKE RIM WITH PLASTIC PUSHRIM 24X 1" SOLID RUBBER TYRE
UPHOLSTERY	BLACK NYLON
WEIGHT	LOAD CAPACITY 100KG GROSS WEIGHT 19KG
PACKING	100 X 34 X 85 (CM)

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za