

CTM HS118 MICRO SCOOTER



Type	Mobility Scooters	Curb Climbing	35mm/1.4"
Model Number	HS-118	Turning Radius	1270mm/50"
Overall Length	1000mm/39.4"	Suspension	N/A
Overall Width	510mm/20"	Brake	Electro-Mechanical
Overall Height	890mm/35"	Seat Type	Padded Foldable
Wheels: Front	195mm/8"	Seat Width	425mm/17"
Wheels: Rear	195mm/8"	Motor Size	250W, 4700 rpm
Weight w/ Batteries	40kg/88lbs	Battery Size	(2) 12V • 12Ah
Max. Speed	6.4kmph/4mph	Battery Weight	9.2kg/20lbs
Weight Capacity	115kg/253lbs	Travel Range	10km/6.2miles
Ground Clearance	30mm/1.2"	Battery Charger	1.8A Off Board
Grade Climbable	8 degrees	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Blue		

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za