

CTM MID-WHEEL DRIVE POWER



Rehab Seat Option

Features

- * Maximum Speed: 6.7kmph / 4.2 mph
- * Capacity: 136kg / 300lbs
- * Full Coli – Over suspension
- * Strong and Responsive Front and Rear Caster Keep you SAFELY on course



STD Seat

NAME	CTM MID-WHEEL DRIVE POWER
PRODUCT CODE	HS2850
OVERALL LENGTH WIDTH	995MM / 39" 610MM / 24"
WEIGHT W/ BATTERIES	80KG / 176LBS
MAX. SPEED	6.7 KMPH / 4.2 MPH
WEIGHT CAPACITY	136KG / 300LBS
GROUND CLEARANCE	60MM / 2.4"
GRADE CLIMBABLE	6 DEGREES
CURB CLIMBING	40MM / 1.6"
TURNING RADIUS	620MM / 24.4"
SEAT WIDTH	420~475~500MM / 17.7"~18.7"~19.6"
ARMREST HEIGHT	225~250~270MM / 8.9"~9.8"~10.6"
DRIVE TRAIN	2-MOTOR MID-WHEEL DRIVE
MOTOR SIZE	420W 4900 RPM
TRAVEL RANGE	20KM / 12.4 MILES
SEAT TYPE	PADDED, BREATHABLE, BALLISTIC NYLON BACK REST

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za