

CTM HS320 MID-LINE SCOOTER



Large Seat Option



Solid Tyre Option

CTM MID-LINE SCOOTER – HS320

| Type | Mobility Scooters | Curb Climbing | 50mm/2" |
|---------------------|--------------------------------------|-----------------|---|
| Model Number | HS-320 | Turning Radius | 1020mm/40.2" |
| Overall Length | 1090mm/42.9" | Suspension | N/A |
| Overall Width | 550mm/21.7" | Brake | Electro-Mechanical |
| Overall Height | 870mm/34.3" | Seat Type | Light Weight Foldable Swivel |
| Wheels: Front | 230mm/9" | Seat Width | 460mm/18" |
| Wheels: Rear | 230mm/9" | Motor Size | 450W 3700 rpm |
| Weight w/ Batteries | 53.8kg/118.7lbs | Battery Size | (2) 12V • 18 Ah |
| Max. Speed | 8kmph/5mph | Battery Weight | 12.7kg/28lbs |
| Weight Capacity | 135kg/300lbs | Travel Range | 135km/8.4 miles |
| Ground Clearance | 40mm/1.6" | Battery Charger | 2A Off Board |
| Grade Climbable | 12 degrees | Electronics | On/Off Key Switch, Battery Level Indicator, Speed Control Knob |
| Colour | Burgundy, Silver, Blue, Turquoise | | |

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za