

# CTM BARIATRIC POWER



Solid Tyre Option

### Features:

- Maximum Speed: 8kmph/5mph
- Capacity: 181 kg / 400lbs
- 500mm / 20" reclining captains chair with seat belt
- Swing-away joystick controller
- ZERO effort single free wheel lever

PRODUCT CODE: HS5600



HS-5600			
Type	Power Chairs	Curb Climbing	50 mm / 2"
Model Number	HS-5600	Turning Radius	790 mm / 31.1"
Overall Length	1135 mm / 44.7"	Supension	N/A
Overall Width	635 mm / 25"	Brake	Electro-Mechanical
Overall Height	1185 mm / 46.7"	Seat Type	Captain Seat, Adjustable Seat Height, & Seat Belt
Wheels: Front	230 mm / 9"	Seat Width	500 mm / 20"
Wheels: Rear	360 mm / 14"	Motor Size	420W, 4600 r.p.m
Weight w/ Batteries	102 kg / 224 lbs	Battery Size	(2) 12V. 50Ah
Max. Speed	8 kmph / 5 mph	Battery Weight	31.3 kg / 69 lbs
Weight Capacity	181 kg / 400 lbs	Travel Range	24 km / 15 Miles
Ground Clearance	85 mm / 3.4"	Battery Charger	5A Off Board
Grade Climbable	8 degree	Electronics	N/A
Color	Burgundy, Silver, Blue		

\* Subject to change without notice.

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
 P.O. Box 51124, Raedene, 2124  
 enquiries@medop.co.za | www.medop.co.za  
 JHB Tel: 011 827 5893/4/5  
 CT Tel: 021 534 8628 | admin@medop.co.za