

# CTM HS570 GT-3 SCOOTER



Solid Tyre Option

## CTM GT-3 SCOOTER – HS570

Type	Mobility Scooters	Curb Climbing	50mm/2"
Model Number	HS-570	Turning Radius	1080mm/42.5"
Overall Length	1240mm/48.8"	Suspension	N/A
Overall Width	610mm/24"	Brake	Electro-Mechanical
Overall Height	965mm/38"	Seat Type	Premium Upholstery Mid-Back w/ Headrest Swivel Seat-Sliding Adjustment
Wheels: Front	260mm/10"	Seat Width	455mm/18"
Wheels: Rear	260mm/10"	Motor Size	500W 3400 rpm
Weight w/ Batteries	84kg/185.2lbs	Battery Size	(2) 12V• 36 Ah
Max. Speed	8kmph/5mph	Battery Weight	24.7kg/55lbs
Weight Capacity	135kg/300lbs	Travel Range	30km/18.6 miles
Ground Clearance	40mm/1.6"	Battery Charger	5A Off Board
Grade Climbable	12 degrees	Electronics	On/Off Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Silver, Blue, Turquoise		

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za