


CTM BARIATRIC SCOOTER



PRODUCT CODE: H5890



Solid Tyre Option

Type	Mobility Scooters	Curb Climbing	70mm/2.8"
Model Number	HS-890	Turning Radius	1730mm/68.1
Overall Length	1530mm/60.2"	Suspension	Front & Rear
Overall Width	685mm/27"	Brake	Electro-Mechanical
Overall Height	1070mm/42.1"	Seat Type	Sliding and Swivel Reclining
Wheels: Front	320mm/13"	Seat Width	480mm/19"
Wheels: Rear	320mm/13"	Motor Size	700W 5400 rpm
Weight w/ Batteries	114kg/251lbs	Battery Size	(2) 12V • 50Ah
Max. Speed	15kmph/9.3mph	Battery Weight	31.3kg/69lbs
Weight Capacity	227kg/500lbs	Travel Range	37.5km/23.4 miles
Ground Clearance	85mm/3.3"	Battery Charger	5A Off Board
Grade Climbable	10 degrees	Electronics	On/Off Key Switch, Battery Level Indicator, Speed Control Knob
Colour	Burgundy, Silver, Blue, Turquoise		

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za