

# DELTA A7 POWER WITH SEAT ELEVATE



<b>NAME</b>	DELTA A7 POWER WITH SEAT RAISER
<b>PRODUCT CODE</b>	A7POW
<b>SIZE</b>	18"
<b>BATTERIES</b>	2X36AH/12 VOLT
<b>MOTORS</b>	25 WATT
<b>CASTORS</b>	8X2" FOAM FILLED GREY
<b>SPEED</b>	6KM/H
<b>RANGE</b>	20KM (TERRAIN AND USER DEPENDENT)
<b>REAR RIMS &amp; TYRES</b>	121/2 X 21/4 ALUMINIUM RIM 121/2 X 21/4 ALUMINIUM RIM
<b>ARMS</b>	HEIGHT & WIDTH ADJUSTABLE ARM
<b>GRADIENT</b>	12 DEGREES
<b>CONTROLS</b>	DYNAMIC SHARK REMD11 POWER MODULE PMB11
<b>SEAT FUNCTION</b>	SEAT ELEVATE MODEL
<b>LOAD CAPACITY</b>	LOAD CAPACITY - 115KG WEIGHT - 85KG
<b>PACKING</b>	71X74X101CM

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za