

GRIFFIN SEMI-ELECTRIC STAND-UP



NAME	SEMI-ELECTRIC STAND-UP
PRODUCT CODE	50003
CONTROLS	SINGLE ACTUATOR, FINGER OPERATED
CASTORS	6X1" GREY CASTORS
REAR RIM	24X1 3/8" PU GREY TYRE ON ALUMINIUM SPOKE RIM
BATTERIES	2 X SEALED 3.8AH/24VOLT
LOAD CAPACITY	LOAD CAPACITY - 110KG WEIGHT - 30 KG
ADJUSTABILITY	BACK HEIGHT, SEAT DEPTH, FOOTREST HEIGHT, CENTRE OF GRAVITY
PACKING SIZES	65X68X94CM 14", 16", 18"

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za