

INVACARE 9000XT FULL RECLINER



Invacare® 9000 XT Recliner Wheelchair Features

Dynamic recline range from 95° to 175°, +/- 5°
Rear axle base moves, providing stable centre-of-gravity in recline position.
Multiple wheel offerings allow variation of seat-to-floor heights
Urethane tread tires mounted on “no flex” wheels offer superior performance
Urethane casters for better ride and lighter weight
Low-maintenance, dual-sealed precision bearings
Heavy-duty inner liners keep seat and back from stretching
Stocked with attendant (standard) or patient-operated wheel locks

NAME	INVACARE 900XT FULL RECLINER
PRODUCT CODE	9000XT
SIZE	14, 16, 18, 20, 22"
FRAME	MILD STEEL BLACK PRISM POWDER COATED
FOOTREST	REMOVABLE
CASTOR	6" / 8" URETHANE CASTORS
REAR RIM	URETHANE REAR TIRES
HEAD REST	10" REMOVABLE HEADREST
UPHOLSTERY	BLACK NYLON
CARRY CAPACITY	LOAD CAPACITY 250lb. WEIGHT - 39lb/18kg
ARMS	FLIP BACK FIXED LENGTH
PACKING	32 X 13 X 37"

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za