

INVACARE ACTION 5



PRODUCT CODE: SJC0110

The Invacare Action5 has been designed for the more active user looking for an affordable chair with good drive performance and the ability to fold easily for handling, storage and transportation purposes. With its two patents and lightweight construction, the new Action5 gives the user a chair they can rely on. Key features include:

- Unique folding mechanism
- Stability of a rigid chair
- Quick step-less adjustability
- Highly configurable

TECHNICAL DATA

Invacare Action ⁵	305 / 330 / 355 / 380 / 405 / 430 / 455 / 480 / 505 / 530 mm	350 / 375 / 400 / 425 / 450 / 475 / 500 mm	0 - 9°	Front 360-535 mm in steps of 25 mm Rear: 360-510 mm in steps 25 mm	Fixed 90° Or stepless 90° -15°/+15°	Fixed 300 / 325 / 350 / 400 / 425 / 450 / 475 / 500 mm Adjustable 360-400 and 410-490 mm	Fixed front 80° Swing away legrest 70° / 80° & 90° Elevating legrest 0-80°
Invacare Action ⁵	Swing away frame: 200-510 mm Fixed frame: 300-530 mm	186-227 mm 186-266 mm 205-280 mm	SW +160 mm To SW+ 250 mm	320 -380 mm	845 -1205 mm	670 -860 mm 701 -1108 mm	
Invacare Action ⁵	7°	1050 -1380 mm	100 & 130 kg	From 12,3 kg	From 8,3 kg		

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za