

INVACARE AZALEA MINOR



PRODUCT CODE: BNC0014

Rea Azalea Minor – designed for petite adults and teenagers

The Rea Azalea Minor is designed for users that require a seat that has a smaller depth and width. Adapted from the Rea Azalea, the Azalea Minor boasts all the advantages of a reliable, tilt-able wheelchair with its special weight-shifting mechanism. This model is ideal for all users that require a solid, stable wheelchair. Now with new frame colours and upholstery. A tension adjustable backrest is available as a new option.



Tension adjustable backrest
A new option adjustable in height and as wide as the seat.



Table tray
It is transparent and the rail offers gripping opportunities for posture solutions.



Shorter Chair
The chassis is 30 mm shorter, the seat depth is 50 mm shorter and the backrest 100 mm lower than on standard models. Here in new colour, electric green.

Technical data

Rea Azalea Minor	340 - 440 mm	380 - 450 mm	400/450 mm	550 - 650 mm	240 - 380 mm	330 - 500 mm
Rea Azalea Minor	0° - 30°	SW + 220 mm (SW + 230 mm)	900 - 1250 mm	900 - 970 mm	32 kg	Max 75 kg Transport weight 14,5 kg
	Rea Azalea Minor complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: www.invacare-reas.com under Rea wheelchairs/Technical downloads					Invacare International Sàrl Route de Cité-Cuest 2 - 1196 Gland Switzerland Tel: +41 (0)22 354 60 10 Fax: +41 (0)22 354 60 11 hqeuropa@invacare.com www.invacare.eu.com
Rea Azalea Minor	Crash-tested ISO 7176-19					

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za