

# INVACARE HEAVY DUTY TRANSPORT CHAIR

**NAME**

INVACARE  
HEAVY DUTY  
TRANSPORT  
CHAIR

**PRODUCT CODE**

TRHD22FR

**SEAT WIDTH**

22"

**SEAT DEPTH**

18"

**BACK HEIGHT**

20"

**SEAT-TO-FLOOR HEIGHT**

19.5"

**OVERALL WIDTH OPEN**

28"

**OVERALL DEPTH W/O RIGGINGS**

36.5"

**LOAD CAPACITY**

400LB

**Features**

- 13 month limited warranty on side frames and cross braces
- 12" rear wheels with flat free tyres
- Removable padded desk arms
- Carbon steel frame with dual cross braces
- Swingaway footrest with aluminium footplates
- Flame retardant black nylon upholstery

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za