

INVACARE STREAM POWER



Maintain independence

The **Stream**, a new power wheelchair from Invacare is a popular model suitable for with reduced mobility. The Invacare **Stream** is ideal for both indoor and outdoor use. To ensure a comfortable sitting position and a relaxing ride, the seat can be adjusted to suit customers unique requirements. A dedicated lighting kit is available as an option, thus extending the possibilities of use for this model.

PRODUCT CODE: 1533501



Light-kit
Available as an option - includes bright headlights with indicators.



Rear wheels
14" pneumatic rear wheels with studded profile.



Swing in/out footrest
The leg rests swing in and out for easy transfer and transportation.

Technical data

Stream	44 - 49 cm ¹	41, 46, 51 cm	42,5 - 50 cm ²	48 cm ²	29 - 36 cm ²	30 - 43 cm ²	up to 22°
Stream	-10° - 30°	64 cm ³	100 cm	103 cm	105 kg	130 kg	50 Ah
Stream	166 cm	5 cm	11%	2 x 180 W	25 km ⁴	Frame colour Gloss aluminium	Upholstery Black nylon
Stream	6 or 9 km/h	<ol style="list-style-type: none"> 1/ seat cushion 43 cm 2/ measured without seat cushion / upholstery 3/ total width is depending on seat adjustment 4/ according to ISO 7176-4:2008 					

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za