

JUMPER POWER WHEELCHAIR



Folding frame Power wheelchair

TECHNICAL SPECIFICATIONS

Product Code	Y205
Material	Aluminium
Sizes	Available in 16", 18", 19"
Wheels	Front Castor 9" foam filled Rear wheel 14" Foam filled
Battery	Lead acid battery 12V 30Amp
Motor	320watt*2
Slope	8 degrees
Armrest	Sloped Height Adjustable Arm
Controller	PG VSI 50A (UK)
Load Capacity	120kg
Colours	Red & Silver
Footrest	Swing Away Footrest
Seat	Detachable quality seat and back cushion, easy to clean
Special:	Advanced height adjustable anti-tip wheels. Highly reliable industry-proven electronics by PG Drives Technology
Max speed	5-8km
Travel mileage	20-25km
Weight of Wheelchair	40kg (excl battery - 20kg)



Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za