

# KUSCHALL COMPACT ATTRACT
















PRODUCT CODE: KUSCO

The küschall Compact attract meets the expectations of users who are looking for a reliable, cross- foldable active wheelchair, offering all the adjustability required for secure seating, maneuverability and easy driving.

A comfortable, practical and efficient wheelchair, the küschall Compact attract comes with swing-away removable leg rests and is fully optimized for transfers. Four matte colours make it an appealing wheelchair for everyone.



## Technical data

 Seat width (SB): 360 - 500 mm in increments of 20 mm	 Seat depth (ST): 360 - 500 mm in increments of 20 mm	 Seat height (SHv, SHh): front: 430 - 520 mm rear: 400 - 490 mm in increments of 10 mm
 Knee-to-heel length (UL): 380 - 460 mm in increments of 10 mm	 Backrest height (RH): 330 - 495 mm in increments of 15 mm	 Backrest angle (RW): 82°/86°/90°/ 94°/98°/102°
 Max. safe slope: 7°	 Reversing width: 920 - 1125 mm	 Total length (GL): 70°: ca. 985 mm 90°: ca. 895 mm w/out legrest: 690 mm
 Total width: SB +180 mm	 Total width folded: ca. 300 mm	 Total weight, unloaded: ca. 11.3 kg
 Total weight w/out rear wheels: ca. 8.7 kg	 Max. user weight: 130 kg	 Crash Test approved Fulfills the requirements of ISO 7176-19



# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za