

# KUSCHALL K-SERIES ATTRACT


















## FEATURES AND BENEFITS

- 1. CONFIGURABLE:** It is available in a wide range of seat sizes and seat heights. Combined with a wide range of options, accessories and colours, the K-Series can be tailored to suit your exact needs and personality.
- 2. ADJUSTABLE:** Adjustments come as standard, so finding your perfect position is guaranteed!  
Possible Adjustments:  
Front seat height                      Rear seat height  
Centre of gravity                      Backrest height and angle
- 3. EASY TRANSPORT:** Quick and easy folding with the perfect open frame design ensures effortless car transfers. Keep independent and on the go!
- 4. MATRX BACKRESTS AND CUSHIONS:** Our range offers additional posture positioning giving you the freedom to enjoy the benefits of life in an active wheelchair.

AN ATTRACTIVE WHEELCHAIR FOR EVERYONE.

## TECHNICAL DATA

 <b>Frame angle</b> 75°/90° V-frame or tapered	 <b>Seat width</b> 320 – 500 mm in increments of 20 mm	 <b>Seat depth</b> 350 – 525 mm in increments of 25 mm	 <b>Seat height</b> front: 450 – 520 mm rear: 380 – 490 mm in increments of 10 mm
 <b>Knee-to-heel length</b> 290 – 480 mm in increments of 10 mm	 <b>Backrest height</b> 270 – 480 mm in increments of 15 mm	 <b>Backrest angle</b> 74°/78°/82°/86°/90°	 <b>Max. safe slope</b> 7°
 <b>Reversing width</b> 920 – 1230 mm	 <b>Total width</b> SW + 170 mm	 <b>Centre of gravity</b> 80 mm/95 mm/110mm/ 125 mm/140 mm	 <b>Total length</b> 75°: approx. 830 - 1055 mm 90°: approx. 760 - 985 mm
 <b>Total weight unloaded</b> approx. 7.8 kg	 <b>Transport weight w/out rear wheels</b> approx. 5.2 kg	 <b>Max. user weight</b> 130 kg	

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za