

ROC HYDRAULIC STAND-UP



Innovate, simple, light weight design to stand-up and sit-down

- * Advanced stability design of safety concern
- * Unique easy control make stand-up and sit-down possible
- * Dual safety protection to prevent accidental improper operation
- * Multiple adjustments to fit user size

NAME	ROC (SMC-S0) HYDRAULIC STAND-UP
PRODUCT CODE	49311
CONTROLS	DUAL GAS ASSISTED LIFT PUMPS
CASTORS	6X1" GREY CASTORS
REAR RIM	24X1 3/8" PNEUMATIC GREY TYRE ON ALUMINIUM RIM
BATTERIES	
LOAD CAPACITY	LOAD CAPACITY - 110KG WEIGHT - 30KG
ADJUSTABILITY	BACK HEIGHT, SEAT DEPTH, FOOTREST HEIGHT, CENTRE OF GRAVITY
PACKING	65X68X94CM
SIZES	14", 16", 18"

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za