

SATURN ALU-LITE KIDDIES



NAME	SATURN ALU-LIGHT
PRODUCT CODE	LW25-77K
SIZE	14" X 14", 14" X 16"
COLOURS	METALLIC BLUE, WET BLACK
FRAME	FOLDING ALUMINIUM POWDER COATED
ARMS	FIXED HEIGHT FLIP BACK FULL ARM
FOOTREST	REMOVABLE
CASTORS	8X1" SOLID CASTOR ON ALUMINIUM FORK
RIMS	24X13/8 MAG RIM WITH PU TYRE 12MM QUICK RELEASE AXLES
UPHOLSTERY	BLACK REINFORCED NYLON
ADJUSTMENT	CENTRE OF GRAVITY, LOW SLUNG HEIGHT ADJUSTABLE BACK REST WITH PUSH HANDLES
WEIGHT	LOAD CAPACITY – 115KG GROSS WEIGHT – 17KG (100X34X85CM)

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za