

SATURN ALU-LITE



NAME	SATURN ALU-LITE
PRODUCT CODE	LW25-77
SIZE	14", 16:, 18" AND 20"
COLOURS	METALLIC BLUE, WET BLACK
FRAME	FOLDING ALUMINIUM POWDER COATED
ARMS	FIXED HEIGHT FLIP BACK DESK ARM
FOOTREST	REMOVABLE
CASTORS	8X1" SOLID CASTOR ON ALUMINIUM FORK
RIMS	24X1 ^{3/8} MAG RIM WITH PU TYRE - 12MM QR AXLE
UPHOLSTERY	BLACK REINFORCED NYLON TAB
ADJUSTMENT	CENTRE OF GRAVITY, LOW SLUNG, HEIGHT ADJUSTABLE BACKREST WITH PUSH HANDLES
OPTIONS	ELEVATING LEG REST, HEIGHT ADJUSTABLE ARMS, 8X2" CASTOR ON STEEL FORK, ALUMINIUM SPOKE RIM, TENSION ADJUSTABLE BACK UPHOLSTERY, AMPUTEE BRACKET, ANTI-TIP WHEELS
WEIGHT	LOAD CAPACITY 115KG GROSS WEIGHT 14 KG
PACKING	100 X 34 X 85 (CM)

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za