

## SATURN COASTER



|                             |  |
|-----------------------------|--|
| <b>NAME</b>                 | SATURN COASTER   |
| <b>PRODUCT CODE</b>         | COA  |
| <b>SIZE</b>                 | 14", 16", 18", 20"   |
| <b>COLOURS</b>              | COBOLT BLUE  |
| <b>FRAME</b>                | REINFORCED POWDER COATED MILD STEEL  |
| <b>ARMS</b>                 | REMOVABLE DESK   |
| <b>FOOTREST</b>             | REMOVABLE, ELEVATING LEG REST UPCHARGE<br>OPTION   |
| <b>CASTORS</b>              | 8x1" SOLID CASTOR ON STEEL FORK  |
| <b>RIMS</b>                 | 24X1 <sup>3/8</sup> MAG RIM WITH 24X1 <sup>3/8</sup> PU TYRE - FIXED<br>AXLE                           |
| <b>UPHOLSTERY</b>           | BLACK REINFORCED NYLON   |
| <b>STRENGTH</b>             | SABS COMPLIANT AND TESTED, CHROME<br>PUSHRIMS, STEEL SIDE PANELS, DOUBLE<br>CROSSBARS ON THE 20" MODEL |
| <b>LOAD CAPACITY/WEIGHT</b> | LOAD CAPACITY 120KG WEIGHT - 20KG  |
| <b>PACKING</b>              | 100 X 35 X 85 (CM)   |

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za