

## SATURN KIDDIES



<b>NAME</b>	SATURN KIDDIES
<b>PRODUCT CODE</b>	KIDDIES
<b>SIZE</b>	12"
<b>COLOURS</b>	WET BLACK
<b>FRAME</b>	FOLDING MILD STEEL POWDER COATED
<b>ARMS</b>	HEIGHT ADJUSTABLE FULL ARM
<b>FOOTREST</b>	REMOVABLE 90° FOOTREST
<b>CASTORS</b>	8X2" PU CASTOR ON STEEL FORK
<b>RIMS</b>	20 X 1.75 STEEL SPOKE RIM WITH 24X1.75 TYRE WITH NO MORE FLAT (NMF), FIXED AXLES
<b>UPHOLSTERY</b>	BLACK REINFORCED NYLON
<b>WEIGHT</b>	LOAD CAPACITY – 50KG GROSS WEIGHT – 16KG
<b>PACKING</b>	100X34X85CM

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za