

SATURN PORTERS



Available with 12" rear rims option

NAME	SATURN PORTER WHEELCHAIR
PRODUCT CODE	POR020 / POR020T
SIZE	20"
COLOURS	BLACK
FRAME	NON FOLDING MILD STEEL POWDER COATED REINFORCING PUSH BAR AND BAR BETWEEN FOOTREST HANGERS
ARMS	FIXED FULL ARM
FOOTREST	FIXED FOOTREST WITH REINFORCING BAR
CASTORS	8X2" PU CASTOR ON 2MM STEEL FORK
RIMS	24X1 ^{3/8} MAG RIM WITH 24X1 ^{3/8} PU TYRE - FIXED AXLE
UPHOLSTERY	BLACK REINFORCED NYLON WITH "Z" FOLD POCKET FOR FILES
STRENGTH	CHROME PUSHRIMS, STEEL SIDE PANELS, RIGID REINFORCED BOX SHAPE DESIGN
WEIGHT	LOAD CAPACITY - 125KG WEIGHT 24KG
PACKING	110 X 80 X97 (CM)
ACCESSORIES	STAINLESS STEEL DRIP STAND "T" ROD
RIMS	12 ½ x 2 ¼ MAG W/Solid tyre

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za