

TOP END CROSSFIRE ALL TERRAIN



The **Invacare® Top End® Crossfire™ All Terrain Wheelchair** is the monster truck of wheelchairs and is designed to leave some serious tracks. Do you like to hike, fish, camp, hunt, kayak, garden or go to the beach? You can do it all with the Crossfire All Terrain Wheelchair. Built of 6061 T6 Aluminium, it is equipped with enormous 25" x 2.75" knobby wheels and huge 8" x 2" front casters for extreme traction. The Crossfire All Terrain Wheelchair is engineered with mega ground clearance, plus forward and side-to-side stability so that you can conquer the toughest terrain.

Features

- 6061 T6 Aluminium frame
- Adjustable Centre-of-gravity
- Camber 0-9 °
- Lock down folding back
- Adjustable back height
- Adjustable tension seat upholstery
- Adjustable tension back upholstery
- 24" or 25" performance Sun® wheels
- Hard anodised handrims
- 24" or 25" x 2.75" knobby tired
- 8" x 2" front casters
- Adjustable height and angle footplate
- Powder coat paint
- Weight capacity 250 lb.



1. 24" or 25" x 2.75" knobby wheels and tires
2. 8" x 2" front casters
3. Engineered for front and side-to-side stability
4. Center-of-gravity adjustment
5. Anti-tip option
6. Optional push handles, wheel locks and anti-tip bars

Medop cc

Medopark, 1 Manchester Road, Wadeville,
 P.O. Box 51124, Raedene, 2124
 enquiries@medop.co.za | www.medop.co.za
 JHB Tel: 011 827 5893/4/5
 CT Tel: 021 534 8628 | admin@medop.co.za