

Manual Wheelchair Low Active

Saturn Coaster
Saturn Hercules
Basik Wheelchair
Basik Aluminium Wheelchair
Basik Hybrid Granny
Swift
Saturn Hercules – XL
Jade Bariatric
Invacare Bariatric IVC 9000 Topaz

SATURN COASTER



NAME	SATURN COASTER
PRODUCT CODE	COA
SIZE	14", 16", 18", 20"
COLOURS	COBOLT BLUE
FRAME	REINFORCED POWDER COATED MILD STEEL
ARMS	REMOVABLE DESK
FOOTREST	REMOVABLE, ELEVATING LEG REST UPCHARGE OPTION
CASTORS	8x1" SOLID CASTOR ON STEEL FORK
RIMS	24X1 ^{3/8} MAG RIM WITH 24X1 ^{3/8} PU TYRE - FIXED AXLE
UPHOLSTERY	BLACK REINFORCED NYLON
STRENGTH	SABS COMPLIANT AND TESTED, CHROME PUSHRIMS, STEEL SIDE PANELS, DOUBLE CROSSBARS ON THE 20" MODEL
LOAD CAPACITY/WEIGHT	LOAD CAPACITY 120KG WEIGHT - 20KG
PACKING	100 X 35 X 85 (CM)

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

SATURN HERCULES



NAME	SATURN HERCULES
PRODUCT CODE	HER
SIZE	14", 16", 18", 20"
COLOURS	COBOLT BLUE
FRAME	REINFORCED POWDER COATED MILD STEEL
ARMS	REMOVABLE DESK OR FULL ARM
FOOTREST	REMOVABLE, ELEVATING LEG REST OPTIONAL EXTRA
CASTORS	8x2" PU CASTOR ON 5MM STEEL FORK
RIMS	24X1 ^{3/8} MAG RIM WITH 24X1.75 PU TYRE - FIXED AXLE
UPHOLSTERY	BLACK REINFORCED NYLON
STRENGTH	SABS COMPLIANT AND TESTED, CHROME PUSHRIMS, STEEL SIDE PANELS, DOUBLE CROSSBARS ON THE 20" MODEL
LOAD CAPACITY/WEIGHT	LOAD CAPACITY 125KG WEIGHT - 21KG
PACKING	100 X 35 X 85 (CM)

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

BASIK WHEELCHAIR



NAME	BASIK WHEELCHAIR
PRODUCT CODE	YK9031
SIZE	16", 18", 20"
COLOURS	BLACK
FRAME	MILD STEEL POWDER COATED
ARMS	FLIP BACK DESK ARMS
FOOTREST	REMOVABLE WITH PLASTIC FOOTPLATE
CASTORS	8X1" SOLID CASTOR
RIMS	24X1" SPOKE RIM WITH PLASTIC PUSHRIM 24X 1" SOLID RUBBER TYRE
UPHOLSTERY	BLACK NYLON
WEIGHT	LOAD CAPACITY 100KG GROSS WEIGHT 19KG
PACKING	100 X 34 X 85 (CM)

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

BASIK ALUMINIUM WHEELCHAIR



YK9060

- Aluminium chair frame
- Flip back Desk Arm
- Tool-free adjustable Swing Away Footrest
- Quick released 24" spoke wheel with solid tire
- 8" front castor

ITEM NO.	SIZE	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to seat Height	Weight Capacity (Kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (Kg)
YK9060	18"	103	64	92	30	46	40	40	48	115	8"	24"	14.5/13

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

BASIK HYBRID GRANNY



- Wheelchair/Transport chair two way use
- Aluminium chair frame
- Flip back desk arm
- Detachable mag wheel with PU tyre
- 8" rear wheel & 8" front castor
- Tool-free adjustable swing away footrest
- Aluminium front fork with 2 adjustable positions

NAME	BASIK HYBRID GRANNY
PRODUCT CODE	Yk9061
SIZE	18"
COLOURS	SILVER
FRAME	ALUMINIUM CHAIR FRAME
ARMS	FLIP BACK DESK ARM
FOOTREST	TOOL-FREE ADJUSTABLE SWING AWAY FOOTREST
CASTORS	24" REAR WHEEL & 8" FRONT CASTOR
RIMS	FIXED MAG WHEEL WITH PU TYRE
UPHOLSTERY	BLACK NYLON
WEIGHT	LOAD CAPACITY 90KG & WEIGHT 16KG
PACKING	100 x 34 x 85 (cm)
USE	LIGHT INDOOR / SMOOTH SURFACE USE

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

SWIFT



BACK UPHOLSTERY



ARM REST



Specification

Item No. SWIFT YK9063-16", 18"

Aluminum frame wheelchair
Powder coating wheelchair frame – silver
Powder coating black cross bars,
24" spoke rear wheel with PU tires
8"x1" PU front castor,
Plastic handrim,
Height/Adj flip up armrest with slope,
Swing in/out & removable footrest,
Plastic front fork,
Plastic footplate,
Plastic side panel
Nylon padded seat
TAB backrest

Carton size: 82x27x93
CMGW/NW: 16.6 / 14.6 KG

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

SATURN HERCULES - XL



NAME	SATURN HERCULES
PRODUCT CODE	HER022
SIZE	22"
COLOURS	COBOLT BLUE
FRAME	REINFORCED POWDER COATED MILD STEEL
ARMS	REMOVABLE DESK OR FULL ARM
FOOTREST	REMOVABLE, ELEVATING LEG REST OPTIONAL EXTRA
CASTORS	8x2" PU CASTOR ON 5MM STEEL FORK
RIMS	24X1 ^{3/8} MAG RIM WITH 24X1.75 PU TYRE - FIXED AXLE
UPHOLSTERY	BLACK REINFORCED NYLON
STRENGTH	SABS COMPLIANT AND TESTED, CHROME PUSHRIMS, STEEL SIDE PANELS, DOUBLE CROSSBARS ON THE 20" MODEL
LOAD CAPACITY/ WEIGHT	LOAD CAPACITY 125KG WEIGHT - 21KG
PACKING	100 X 35 X 85 (CM)

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

JADE BARIATRIC



Specifications:

PRODUCT CODE: JADE

Size	22", 24", 26"
Seat Depth	18"
Frame	Reinforced Powder Coated Mild Steel
Footrest	Removable
Upholstery	Black Reinforced Nylon
Load Capacity	145kg – 185kg
Optional	Elevating leg rest

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

INVACARE BARIATRIC IVC 9000 TOPAZ



HDPG option 1, 000 Fb
Product weight (9TPZA)
Limited warranty 3 years

Specifications:

PRODUCT CODE: 9TPZ

Seat Sizes:

Widths 20", 22", 24", 26" 28", 30"

Depths 18", 20"

Seat-to-floor heights 17.5", 19.5"

15.5" with different wheel and caster packages

Back heights – adjustable 16"-18"

Front riggings

Swingaway footrest

Swingaway padded elevating legrest

Arm styles

Full-length fixed height, padded, removable

Desk-length fixed height, padded, removable

Full-length adjustable height, padded, removable

Desk-length adjustable height, padded, removable

(All arm styles available conventional)

Weight capacity 700 lb

HDPG option 1,000 lb

Product weight (9TPZ/WD2818/28) 85lb

Limited warranty 3 years on frame and cross braces

The wheelchair is designed specifically to suit the unique needs of bariatric patients. The Topaz combines the strength and durability of a triple chrome-plated, carbon steel frame with the comfort of heavy-duty front riggings and arms for proper positioning. With a variety of seat widths and depths, frame reinforcements that include gusseted sides, caster journals and dual heavy-duty cross braces, this chair provides not only ease in manoeuvrability but stability that will be trusted for years to come.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za