

BASIK HYBRID GRANNY



- Wheelchair/Transport chair two way use
- Aluminium chair frame
- Flip back desk arm
- Detachable mag wheel with PU tyre
- 8" rear wheel & 8" front castor
- Tool-free adjustable swing away footrest
- Aluminium front fork with 2 adjustable positions

| | |
|---------------------|--|
| NAME | BASIK HYBRID GRANNY |
| PRODUCT CODE | Yk9061 |
| SIZE | 18" |
| COLOURS | SILVER |
| FRAME | ALUMINIUM CHAIR FRAME |
| ARMS | FLIP BACK DESK ARM |
| FOOTREST | TOOL-FREE ADJUSTABLE SWING AWAY FOOTREST |
| CASTORS | 24" REAR WHEEL & 8" FRONT CASTOR |
| RIMS | FIXED MAG WHEEL WITH PU TYRE |
| UPHOLSTERY | BLACK NYLON |
| WEIGHT | LOAD CAPACITY 90KG & WEIGHT 16KG |
| PACKING | 100 x 34 x 85 (cm) |
| USE | LIGHT INDOOR / SMOOTH SURFACE USE |

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za