

INVACARE HEAVY DUTY TRANSPORT CHAIR



NAME

INVACARE
HEAVY DUTY
TRANSPORT
CHAIR

PRODUCT CODE

TRHD22FR

SEAT WIDTH

22"

SEAT DEPTH

18"

BACK HEIGHT

20"

SEAT-TO-FLOOR HEIGHT

19.5"

OVERALL WIDTH OPEN

28"

OVERALL DEPTH W/O RIGGINGS

36.5"

LOAD CAPACITY

400LB

Features

- 13 month limited warranty on side frames and cross braces
- 12" rear wheels with flat free tyres
- Removable padded desk arms
- Carbon steel frame with dual cross braces
- Swingaway footrest with aluminium footplates
- Flame retardant black nylon upholstery

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za