

# WHEELCHAIR ACCESSORIES AND OPTIONS

## SINGLE PIECE RAMPS



PRODUCT CODE: RAM

Aluminium Frame  
 Sizes - 36", 48", 72"  
 Single piece foldable with hinge and handle  
 250kg Load Capacity  
 Width - 34"

## TELESCOPIC RAMPS



PRODUCT CODE: RAM020

15cm Internal width  
 25cm External width  
 3 Section Aluminium Frame  
 Telescopic 110cm to 280cm  
 Gross Weight per set - 25kg  
 5cm Lip Height  
 Rough Grip along the surface and on the lip  
 Load Capacity - 200kg

## CHIN CONTROL WITH SEAT



## TRANSFER BOARD

PRODUCT CODE: TRB002



## EZ ENABLER TRAY

PRODUCT CODE: 10858



ABS table for active wheelchairs  
 Aluminium clamp with e-clip x 2  
 With bag  
 Packing: Box size: 43 x 25.2 x 13cm

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
 P.O. Box 51124, Raedene, 2124  
 enquiries@medop.co.za | www.medop.co.za  
 JHB Tel: 011 827 5893/4/5  
 CT Tel: 021 534 8628 | admin@medop.co.za