

SATURN HERCULES - XL



NAME	SATURN HERCULES
PRODUCT CODE	HER022
SIZE	22"
COLOURS	COBOLT BLUE
FRAME	REINFORCED POWDER COATED MILD STEEL
ARMS	REMOVABLE DESK OR FULL ARM
FOOTREST	REMOVABLE, ELEVATING LEG REST OPTIONAL EXTRA
CASTORS	8x2" PU CASTOR ON 5MM STEEL FORK
RIMS	24X1 ^{3/8} MAG RIM WITH 24X1.75 PU TYRE - FIXED AXLE
UPHOLSTERY	BLACK REINFORCED NYLON
STRENGTH	SABS COMPLIANT AND TESTED, CHROME PUSHRIMS, STEEL SIDE PANELS, DOUBLE CROSSBARS ON THE 20" MODEL
LOAD CAPACITY/ WEIGHT	LOAD CAPACITY 125KG WEIGHT - 21KG
PACKING	100 X 35 X 85 (CM)

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za