

# SATURN HERCULES



<b>NAME</b>	SATURN HERCULES
<b>PRODUCT CODE</b>	HER
<b>SIZE</b>	14", 16", 18", 20"
<b>COLOURS</b>	COBOLT BLUE
<b>FRAME</b>	REINFORCED POWDER COATED MILD STEEL
<b>ARMS</b>	REMOVABLE DESK OR FULL ARM
<b>FOOTREST</b>	REMOVABLE, ELEVATING LEG REST OPTIONAL EXTRA
<b>CASTORS</b>	8x2" PU CASTOR ON 5MM STEEL FORK
<b>RIMS</b>	24X1 <sup>3/8</sup> MAG RIM WITH 24X1.75 PU TYRE - FIXED AXLE
<b>UPHOLSTERY</b>	BLACK REINFORCED NYLON
<b>STRENGTH</b>	SABS COMPLIANT AND TESTED, CHROME PUSHRIMS, STEEL SIDE PANELS, DOUBLE CROSSBARS ON THE 20" MODEL
<b>LOAD CAPACITY/WEIGHT</b>	LOAD CAPACITY 125KG WEIGHT - 21KG
<b>PACKING</b>	100 X 35 X 85 (CM)

**Medop** cc

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za