

BASIK KIDDIES



NAME	BASIK KIDDIES
PRODUCT CODE	YK9031K
SIZE	12", 14"
COLOURS	MULTICOLOURED FRAME
FRAME	FOLDING MILD STEEL POWDER COATED
ARMS	DESK, FLIP BACK
FOOTREST	REMOVABLE
CASTORS	6" CASTOR
RIMS	20" REAR WHEEL
UPHOLSTERY	BLACK NYLON
WEIGHT	LOAD CAPACITY - 60KG GROSS WEIGHT - 18.5KG
PACKING	80X34X85CM

SATURN KIDDIES



NAME	SATURN KIDDIES
PRODUCT CODE	KIDDIES
SIZE	12"
COLOURS	WET BLACK
FRAME	FOLDING MILD STEEL POWDER COATED
ARMS	HEIGHT ADJUSTABLE FULL ARM
FOOTREST	REMOVABLE 90° FOOTREST
CASTORS	8X2" PU CASTOR ON STEEL FORK
RIMS	20 X 1.75 STEEL SPOKE RIM WITH 24X1.75 TYRE WITH NO MORE FLAT (NMF), FIXED AXLES
UPHOLSTERY	BLACK REINFORCED NYLON
WEIGHT	LOAD CAPACITY – 50KG GROSS WEIGHT – 16KG
PACKING	100X34X85CM

SATURN ALU-LITE KIDDIES



NAME	SATURN ALU-LIGHT
PRODUCT CODE	LW25-77K
SIZE	14" X 14", 14" X 16"
COLOURS	METALLIC BLUE, WET BLACK
FRAME	FOLDING ALUMINIUM POWDER COATED
ARMS	FIXED HEIGHT FLIP BACK FULL ARM
FOOTREST	REMOVABLE
CASTORS	8X1" SOLID CASTOR ON ALUMINIUM FORK
RIMS	24X13/8 MAG RIM WITH PU TYRE 12MM QUICK RELEASE AXLES
UPHOLSTERY	BLACK REINFORCED NYLON
ADJUSTMENT	CENTRE OF GRAVITY, LOW SLUNG HEIGHT ADJUSTABLE BACK REST WITH PUSH HANDLES
WEIGHT	LOAD CAPACITY – 115KG GROSS WEIGHT – 17KG (100X34X85CM)

INVACARE ACTION 3 JUNIOR



NAME	INVACARE ACTION 3 JUNIOR
PRODUCT CODE	CNB0030
SIZE	CHILD 8"-1" JUNIOR 12"-15"
COLOURS	OPTIONAL
FRAME	FOLDING ALUMINIUM FRAME
ARMS	REMOVABLE & HEIGHT ADJUSTABLE "T" ARM
FOOTREST	REMOVABLE
CASTORS	6" SOLID CASTER
RIMS	20" ALUMINIUM SPOKE RIM WITH PU TYRE QUICK RELEASE AXLES
UPHOLSTERY	BLACK REINFORCED NYLON
ADJUSTMENT	CENTRE OF GRAVITY LOW SLUNG
WEIGHT	LOAD CAPACITY – 80KG GROSS WEIGHT – 15KG

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

ACTION 3 JUNIOR EVOLUTIF

A wheelchair that grows with the child

Invacare **Action³ Junior** is a lightweight, foldable paediatric wheelchair designed for children aged between 3 and 15 years. **Action³ Junior** has been developed to match the individual needs of the child and can grow as they grow. Because a child's needs will change the **Action³ Junior** offers a large range of options to accompany the child in their development and to provide the necessary clinical support. The **Action³ Junior** comes in attractive, modern styles, designed to appeal to kids.



Seat-width adjustment

The seat width of the chassis can be adjusted with an innovative telescopic crossbar for four fixed-seat positions. (Child: 20,5 to 28 cm, Young: 30,5 to 38 cm.)

Upholstery adjustment

As the chassis is made wider, so the upholstery on the back and on the seat of the chair can be made wider with Velcro flaps.

Table tray

In addition to the standard slide-on table tray, we offer a swing-away tray which can easily be folded away.

Adjustable seat height and balance

For improved driving performance and different height options, the rear wheel can be mounted in 20 different positions.

Technical data

Invacare Action³ Junior Child	20,5 / 23 / 25,5 / 28 cm	20 - 25 / 25 - 30 / 30 - 35 cm	41 / 43,5 / 46 cm	37 / 40 / 43 cm ¹	15 - 25 ² / 18 - 28 cm ³ 13 - 27 ⁴ / 22 - 32 ⁵ / 30 - 45 cm ⁶	
Invacare Action³ Junior Young	30,5 / 33 / 35,5 / 38 cm	27,5 - 32,5 / 32,5 - 37,5 cm	41 / 43,5 / 46 cm	37 / 40 / 43 cm ¹	15 - 25 ² / 18 - 28 cm ³ 13 - 27 ⁴ / 22 - 32 ⁵ / 30 - 45 cm ⁶	
Invacare Action³ Junior Child	3° / ±1°	7° / 0° - 30° ⁷	SW + 19 / + 17,5 cm ⁸	87 - 113 cm	67 - 97 cm	32,5 cm
Invacare Action³ Junior Young	3° / ±1°	7° / 0° - 30° ⁷	SW + 19 / + 17,5 cm ⁸	87 - 113 cm	67 - 97 cm	35 cm
Invacare Action³ Junior Child	12,3 kg	80 kg	7,5 kg			
Invacare Action³ Junior Young	15 kg	80 kg	9 kg			

1/ Look in the EPF for possible adjustment of upholstery and pushhandles

2/ Desk

3/ Junior

4/ 90 degree hanger

5/ 80 degree short hanger

6/ 80 degree long hanger

7/ Recliner

8/ Transit Version

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

MATRIX MINI



Transforming kids seating

The *Invacare Matrix Mini Back* is contoured to fit children of all shapes and sizes, combining exceptional support in an appealing low profile package. The *Mini Back* features Mini Set hardware, a specially designed miniature version of the popular Easy Set hardware. Every element of this back support has been carefully crafted and constructed, keeping pediatric users at the forefront at every stage.



Technical data

Mini	Height 18 cm, 25 cm, 30 cm & 36 cm	Width 25 cm - 36 cm	Depth* Adjustment 10 cm	Angle Adjustment 40° range +/- 20°	Product Weight 0,8 kg and up
Mini	Weight Capacity 68 kg	Limited Warranty 2 Years	* Limitation of depth/angle adjustment if used on chairs narrower than model width		

Hardware

This "Mini Set" hardware installs and adjusts in minutes. There is a small attachment point for minimal interference on the back canes and it can be easily adjusted for growth. The adjustable mounting pins allow it to fit to chairs with the same or 1" wider back. It also has a quick release mechanism for convenient portability.



Lightweight back shell

The durable aluminum shell has been designed for easy attachment of the optional headrest.



Water-resistant, breathable cover

Cover and foam overlap the shell for comfort and protection.



Lumbo Sacral Support

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za

TIAN PAEDS STANDING FRAME

PRODUCT CODE: TC402 & TC404



Code: TC402



Code: TC404

SPECIFICATION

Item	3years old	7 years old	Notes
Foot pad	Adjust to the highest (1220mm)	Adjust to the lowest (1020mm)	<i>The distance to the foot-pad</i>
Knee pad	200	415	
Belly pad	300	700(inversion)	
Chest pad	600	1000	

Attentions and applications

1. This product is suitable for children with cerebral palsy aged 3~7 years
2. User height is shorter than 140cm & weight less than 70kg.
3. Indoor Use.: Net weight:25kg.

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za